

THE MOST COMFORTABLE RECLINER IN THE WORLD!



| COLLECTION            | ULTRA <sup>1,2,3</sup>       | ULTRA <sup>1,2,3</sup>  | ULTRA <sup>1,2,3</sup>     | RADIANCE <sup>1,2,3</sup>    | RADIANCE <sup>1,2,3</sup>    | RADIANCE <sup>1,2,3</sup> | RADIANCE <sup>1,2,3</sup>    | ATLAS PLUS 2 <sup>1,2</sup> | TRANQUIL 2 <sup>1</sup>      | TRANQUIL 2 <sup>1</sup> | TRANQUIL 2 <sup>1</sup>       | ELEGANCE 2 <sup>1</sup> | ELEGANCE 2 <sup>1</sup>       | METRO 2 <sup>1</sup>         | URBANA 2 <sup>1</sup>        |
|-----------------------|------------------------------|-------------------------|----------------------------|------------------------------|------------------------------|---------------------------|------------------------------|-----------------------------|------------------------------|-------------------------|-------------------------------|-------------------------|-------------------------------|------------------------------|------------------------------|
| MODEL                 | PLR-4955S                    | PLR-4955M               | PLR-4955L                  | PLR-3955PW                   | PLR-3955S                    | PLR-3955M                 | PLR-3955LT                   | PLR-2985M                   | PLR-935S                     | PLR-935M                | PLR-935LT                     | PLR-975M                | PLR-975L                      | PLR-925M                     | PLR-965M                     |
| SIZE                  | Small                        | Medium                  | Large                      | Petite Wide                  | Small                        | Medium                    | Large/Tall                   | Medium                      | Small                        | Medium                  | Large/Tall                    | Medium                  | Large                         | Medium                       | Medium                       |
| WEIGHT CAPACITY       | 400 lbs.                     | 400 lbs.                | 400 lbs.                   | 400 lbs.                     | 400 lbs.                     | 400 lbs.                  | 400 lbs.                     | 400 lbs.                    | 400 lbs.                     | 400 lbs.                | 400 lbs.                      | 400 lbs.                | 400 lbs.                      | 400 lbs.                     | 400 lbs.                     |
| POSITION <sup>1</sup> | True-Infinite 5'4" and below | True-Infinite 5'5" - 6' | True-Infinite 6' and above | True-Infinite 5'4" and below | True-Infinite 5'4" and below | True-Infinite 5'5" - 6'1" | True-Infinite 6'2" and above | True-Infinite 5'4"-5'9"     | True-Infinite 5'3" and below | True-Infinite 5'4"-5'9" | True-Infinite 5'10" and above | True-Infinite 5'4"-5'9" | True-Infinite 5'10" and above | True-Infinite 5'3" and below | True-Infinite 5'3" and below |
| OVERALL WIDTH         | 37"                          | 39"                     | 39"                        | 41"                          | 37"                          | 39"                       | 40"                          | 39"                         | 33.75"                       | 37.5"                   | 38"                           | 35.5"                   | 37.5"                         | 33"                          | 33"                          |
| CHAIR WEIGHT          | 158 lbs.                     | 158 lbs.                | 176 lbs.                   | 161 lbs.                     | 156 lbs.                     | 158 lbs.                  | 178 lbs.                     | 164 lbs.                    | 137 lbs.                     | 145 lbs.                | 155 lbs.                      | 152 lbs.                | 157 lbs.                      | 142 lbs.                     | 145 lbs.                     |
| CHAIR HEIGHT          | 41"                          | 43"                     | 47"                        | 40"                          | 41"                          | 43"                       | 49"                          | 42"                         | 40"                          | 42"                     | 47"                           | 42"                     | 43"                           | 41"                          | 42"                          |
| SEATING TYPE          | Total Comfort                | Total Comfort           | Total Comfort              | Total Comfort                | Total Comfort                | Total Comfort             | Total Comfort                | Total Comfort               | Total Comfort                | Total Comfort           | Total Comfort                 | Total Comfort           | Total Comfort                 | Total Comfort                | Total Comfort                |
| SEAT WIDTH            | 20"                          | 22.5"                   | 23"                        | 24"                          | 20"                          | 22.5"                     | 23.5"                        | 22"                         | 20"                          | 22"                     | 23.5"                         | 21"                     | 22"                           | 21"                          | 21"                          |
| SEAT DEPTH            | 18"                          | 21"                     | 22"                        | 18"                          | 18"                          | 21"                       | 22"                          | 21"                         | 18.5"                        | 20.5"                   | 22.5"                         | 21"                     | 22"                           | 21"                          | 20"                          |
| SEAT TO FLOOR HEIGHT  | 18"                          | 19.5"                   | 21.5"                      | 19.5"                        | 18"                          | 19.5"                     | 21.5"                        | 19.5"                       | 18"                          | 19.5"                   | 21.5"                         | 20"                     | 20.5"                         | 19"                          | 20"                          |
| TOP OF BACK TO SEAT   | 25"                          | 28"                     | 30"                        | 23"                          | 25"                          | 28"                       | 31"                          | 26"                         | 25"                          | 27"                     | 30"                           | 27"                     | 29"                           | 26.5"                        | 25"                          |
| DISTANCE FROM WALL    | 15"                          | 15"                     | 21.5"                      | 14"                          | 15"                          | 15"                       | 21.5"                        | 16"                         | 13"                          | 16"                     | 16"                           | 14"                     | 16"                           | 14"                          | 14"                          |
| FULLY PADDED CHAISE   | Yes                          | Yes                     | Yes                        | Yes                          | Yes                          | Yes                       | Yes                          | Yes                         | Yes                          | Yes                     | Yes                           | Yes                     | Yes                           | Yes                          | Yes                          |
| FOOTREST EXTENSION    | Standard                     | Standard                | Standard                   | Standard                     | Standard                     | Standard                  | Standard                     | Standard                    | Standard                     | Standard                | Standard                      | Standard                | Standard                      | Standard                     | Standard                     |
| HEAT                  | Standard                     | Standard                | Standard                   | Standard                     | Standard                     | Standard                  | Standard                     | NA                          | NA                           | NA                      | NA                            | NA                      | NA                            | NA                           | NA                           |
| WARRANTY              | 7-Yr. Pro-rated              | 7-Yr. Pro-rated         | 7-Yr. Pro-rated            | 7-Yr. Pro-rated              | 7-Yr. Pro-rated              | 7-Yr. Pro-rated           | 7-Yr. Pro-rated              | 7-Yr. Pro-rated             | 7-Yr. Pro-rated              | 7-Yr. Pro-rated         | 7-Yr. Pro-rated               | 7-Yr. Pro-rated         | 7-Yr. Pro-rated               | 7-Yr. Pro-rated              | 7-Yr. Pro-rated              |

| STANDARD FABRIC SELECTION GUIDE | CAPRICCIO  | SORRENTO | TRENTO <sup>4</sup> | CANYON |        |          |          | BADLANDS | ASTRO    |          |          | BADLANDS | SAVILLE | STONEWASH |
|---------------------------------|------------|----------|---------------------|--------|--------|----------|----------|----------|----------|----------|----------|----------|---------|-----------|
|                                 |            |          |                     |        |        |          |          |          |          |          |          |          |         |           |
|                                 | SLATE      | COFFEE   | BROWN               | OCEAN  | OCEAN  | MUSHROOM | MUSHROOM | GREY     | GREY     | GREY     | GREY     | MUSHROOM | BROWN   | GRANITE   |
|                                 |            |          |                     | SILT   | SILT   | WALNUT   | WALNUT   | BROWN    | BROWN    | BROWN    | BROWN    | WALNUT   |         | GUNMETAL  |
|                                 | CAPPUCCINO |          | CHARCOAL            | STEEL  | STEEL  | STEEL    | STEEL    | MUSHROOM | MUSHROOM | MUSHROOM | MUSHROOM | STEEL    |         |           |
|                                 |            |          |                     |        | WALNUT |          |          |          |          |          |          |          |         |           |

STANDARD FEATURES:



POWER HEADREST



POWER LUMBAR



USB REMOTE



LITHIUM BATTERY BACKUP



FOOTREST EXTENSION

POSITION GUIDE:



RELAXED



RECLINED



UPWARD LIFT



TRUE INFINITE



= Indicates shortest available lead time.

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.  
2. The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.  
3. Unit has three different heating pads in the shoulders, lumbar and seat.  
4. Available for the PLR-4955M and PLR-4955L only.

\*Not all units achieve all positions. Please refer to specific product.  
\*\*Height ranges should be used as a guide only due to varying factors of each individual user.  
Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.