


# Celebrity® DX

Weight capacity:	159 kg
Maximum speed:	Up to 6 km/h
Ground clearance:	9 cm
Turning radius:	150 cm
Overall length:	119 cm
Overall width:	61 cm
Tyres:	
Front:	9 cm x 25 cm
Rear:	10 cm x 25 cm
Anti-tips:	Rear
Suspension:	Front & rear
Transport weight	
Total w/o batteries:	68.5 kg
Heaviest piece when disassembled:	30 kg
Std. Seating:	Deluxe seat with headrest and manual recline
Battery charger:	Off-board, 5 amp
Per-charge range:	Up to 45 km
Battery requirements:	(2) 12V deep cycle, 40 AH standard
Battery size:	40 AH: 14.5 kg each
Warranty:	2-year limited warranty

Colours:



## FEATURES

- State-of-the-art rear suspension
- Top speed of up to 6 km/h
- Maximum weight capacity of 159 kg
- Ultra-simple, auto-latching lockup for quick disassembly and assembly
- Easily accessible tie down points (for transport of unoccupied scooter)
- Front basket and rearview mirror standard
- Tight turning radius of 150 cm
- Includes standard angle-adjustable light
- Front frame-mounted seat post for optimal stability
- High visibility reflectors
- Microprocessor-based controller
- Tiller rotation lock permits tiller to be used as a handle during transport

## OPTIONS

- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Walker holder



**Delta tiller with wraparound handles**

The information contained herein is correct at the time of publication; due to the manufacturers commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

# CELEBRITY® DX

SUPERIOR COMFORT,  
PERFORMANCE, AND STYLE!

HMV-Nr. 18.51.05.1011

## CELEBRITY® DX

The **Celebrity® DX** delivers superb performance and ease of use with sleep, sporty styling for an overall package that impresses with elegant simplicity. Loaded with standard features like a deluxe seat with headrest and manual recline, and lower light for pathway illumination, the Celebrity DX is equipped to excite.



**Pride**®  
Mobility Products BV  
— LIVE YOUR BEST™ —



+31 251 377 577 (Tel) • +31 251 377 580 (Fax) • [www.pridemobility.com](http://www.pridemobility.com)