

When it comes to your busy, on-the-go lifestyle, you need a scooter that is up to the task. Pride's Go-Go[®] Travel Mobility series¹ offers a wide range of mobility scooters to suit your needs. Enjoy total convenience and exceptional features with Pride's travel scooter lineup!

Most models feature...

- Quick and easy feather-touch disassembly
- Delta tiller with port to conveniently charge your scooter
- Standard front basket



FDA Class II Medical Devices¹



SCAN THE QR CODE
TO LEARN MORE ABOUT
TRAVEL MOBILITY!



FOLLOW US ON SOCIAL MEDIA

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]



GO-GO® ES 2¹



GO-GO® ENDURANCE Li¹



GO-GO® ULTRA X¹



GO-GO ELITE TRAVELLER®¹



GO-GO ELITE TRAVELLER® PLUS¹



GO-GO® LX WITH CTS¹



GO-GO® SPORT¹



GO-GO® FOLDING SCOOTER¹

PART #	COLORS	WEIGHT CAPACITY	TURNING RADIUS ³	MAX RANGE PER CHARGE ^{2,6}	TOP SPEED ⁶	CODING
S81 (3W)	Black	250 lbs.	31"	Up to 6.24 miles	Up to 4 mph	K0800 Group 1 Standard
S54LXLIT	Garnet Red (matte), Euro Grey	325 lbs.	48.25"	8Ah: Up to 7.2 miles ⁴ 16Ah MAXX: Up to 14.4 miles ⁸	Up to 4.2 mph	K0899 Power Mobility Device, Not Coded by DME PDAC or does not meet criteria
S39 (3W) S49 (4W)	Black	300 lbs.	36.5" (3W) 52.25" (4W)	3W: Up to 8 miles (200 lbs.) 3W: Up to 6.2 miles (300 lbs.) 4W: Up to 8.9 miles (200 lbs.) 4W: Up to 6.6 miles (300 lbs.)	3W: Up to 4.3 mph; 4W: Up to 4.1 mph	K0899 Power Mobility Device, Not Coded by DME PDAC or does not meet criteria
SC40E (3W) SC44E (4W)	Red, Blue	300 lbs.	33" (3W) 44" (4W)	3W: Up to 7.25 miles (12A) 3W: Up to 10.5 miles (18A) 4W: Up to 6 miles (12A) 4W: Up to 9.75 miles (18A)	Up to 4 mph	K0800 Group 1 Standard
SC53 (3W) SC54 (4W)	Red, Blue	300 lbs.	34.5" (3W) 45.75" (4W)	3W: Up to 10.25 miles 4W: Up to 14.5 miles	Up to 4.5 mph	K0806 Group 2 Standard
S50LX (3W) S54LX (4W)	Red, Blue, Orion Yellow, Onyx Black, Pearl Pink, and Sunburst Orange ⁷	300 lbs.	33.5" (3W) 44.75" (4W)	3W: Up to 6 miles 4W: Up to 5.5 miles	Up to 4.4 mph	Power Mobility Device, Not Coded by DME PDAC or does not meet criteria
S73 (3W) S74 (4W)	Red, Blue	325 lbs.	36" (3W) 50.5" (4W)	3W: Up to 12.75 miles 4W: Up to 12 miles	Up to 4.7 mph	K0801 Group 1 Heavy Duty
SLA: S19WH1003 Lithium-Ion: S19WH1004	White	250 lbs.	32"	SLA: Up to 9.3 miles Lithium-ion: Up to 13 miles ⁴	Up to 4 mph	K0899 Power Mobility Device, Not Coded by DME PDAC or does not meet criteria

GO-GO® FOLDING SCOOTER

FRAME: 2-year **CONTROLLER:** 18 months
ELECTRONICS & CHARGER: 1-year

ALL OTHER GO-GO® MODELS

FRAME: Lifetime limited **ELECTRONICS:** 2-year limited
DRIVETRAIN: 2-year limited
GO-GO® ENDURANCE Li BATTERY: 2-year limited

WARRANTY

ACCESSORIES⁵

- Cane clip attachment
- Cane/crutch holder
- Cup holder
- Forearm crutch holder
- 50" Lap belt
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger

1) These are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
2) Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
3) Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
4) Lithium-ion battery is subject to federal regulations and airline policy. Consult airline prior to flight.

5) Not all accessories are available for all models.
6) Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This can be subject to a variance of (+ or -) 10%.
7) Red and blue interchangeable shroud panels on the Go-Go® LX with CTS are standard. Orion Yellow, Onyx Black, Pearl Pink, and Sunburst Orange shroud panels are optional.
8) 16AH MAXX battery pack is not airline compliant.