

Model	Go Chair®
Part Number	GoChair
Weight Capacity	300 lbs.
Maximum Speed ¹	Up to 3.8 mph
Ground Clearance ³	1.75" at motor
Turning Radius ³	25.5"
Overall Length ³	28.5" (33" with footplate)
Overall Width ³	25.5"
Front Wheels	6" solid casters
Drive Wheels	9" solid
Rear Casters	2" solid anti-tips
Suspension	Limited
Range Per Charge ^{1,2}	Up to 11.2 miles
Heaviest Piece When Disassembled ³	36 lbs. (rear section)
Standard Seating	Type: Compact, foldable Material: Black vinyl Dimensions: Width: 18", Depth: 17", Weight: 30 lbs.; Width: 20", Depth: 18", Weight: 32 lbs.;
Dual Braking System	Regenerative and electromechanical
Standard Electronics	40 amp, PG GC3
Battery Requirements ⁴	20 Ah battery pack Weight: 31.4 lbs.
Battery Charger	Off-board, 2 amp
Warranty	Frame: 5-year limited Drivetrain: 13-month limited Electronics: 13-month limited Batteries: 13-month limited

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.

PRIDE
MOBILITY
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |      

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

032024

GOCHAIR®

**GET
OUT
THERE!**



FDA Class II
Medical Device¹

LIVE YOUR BEST®

pridemobility.com |      

GET OUT THERE

WITH THE #1 PORTABLE POWER CHAIR!

Convenient storage pocket

Five-point swivel seat with latch release

One-hand feather-touch disassembly

FDA Class II Medical Device¹



GC3 controller with quick disconnect

Dual swing-away, under seat storage bins

Large footplate for added comfort

The Go Chair® is #1 for portability! With convenient feather-touch disassembly, the Go Chair allows you to experience lightweight travel and independence on the go. Enjoy dual swing-away storage bins, an array of contemporary colors, and a compact turning radius for outstanding maneuverability.

FEATURES

- 18" x 17" or 20" x 18" standard 5-point swivel seat with latch release
- Compact turning radius
- Dual swing-away, under seat storage bins
- Height and width adjustable armrests
- One hand feather-touch disassembly
- 60" lap belt
- The heaviest piece is only 36 lbs.

OPTIONS

- Medium back memory seat options available in 18" x 18" or 20" x 18"
- Inline right and left swing-away joystick

AVAILABLE COLORS



ACCESSORIES*

- Cane/crutch holder (3" single, 3" double)
- Cup holder
- Forearm crutch holder
- Cell phone holder
- Lap belt (50", 70")
- Oxygen tank holder
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Additional battery pack

MILITARY PATCHES



* Only one rear accessory can be used at a time.



New GC3 controller with quick disconnect.



Store items close at hand with convenient under seat storage.



Personalize your Go Chair® with one or more optional accessories.