

Part Number	SC443EMS
Weight Capacity	300 lbs.
Maximum Speed²	Up to 4.23 mph
Ground Clearance⁴	1.29" at motor 3.5" under deck
Turning Radius⁴	39.2"
Overall Length⁴	41.5"
Overall Width^{4,6}	19.45"
Seat-to-Deck	16.38 – 17.83"
Seat-to-Ground	20.75 - 22.25
Tires	Front: 2"x 7" Rear: 2.5"x 8"
Range^{2,3}	12 Ah: Up to 7.1 miles at 300 lbs. Up to 8.9 miles at 220 lbs. 18 Ah: Up to 11.3 miles at 300 lbs. Up to 14.1 miles at 220 lbs.
Product Weight	83.6 lbs. without batteries 103 lbs. with 12 Ah batteries 111 lbs. with 18 Ah batteries
Heaviest Piece⁴	36.2 lbs. (front section)
Standard Seating	Type: Foldable, mesh back Material: Black Vinyl Dimensions: 17" (W) x 17" (D) Weight: 23 lbs.
Motor Output	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative & electromechanical
Batteries	(Two) 12-volt, sealed lead-acid, non-spillable deep-cycle
Battery Weight⁵	12 Ah: 20 lbs. (standard) 18 Ah: 27.4 lbs.
Battery Charger	3.5 amp off-board charger
Warranty	Lifetime limited on frame; 2-year limited on electronics; 2-year limited on drivetrain 6 months on battery

HCPCS Code: K0899

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 4 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 5 Battery weight may vary based on manufacturer.
- 6 Dimensions listed are for scooter only and do not include seats or arms. Width may vary based on seat options and accessories.



Featuring a new
deck design with
more leg room!



FDA Class II
Medical Device

PRIDE[®]
MOBILITY

pridemobility.com | 800.800.8586

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

GO GO[®]
ELITE TRAVELLER GEN 3
PLATINUM



FDA Class II
Medical Device



pridemobility.com

GO GO[®]

ELITE TRAVELLER GEN 3

PLATINUM

FDA Class II Medical Device

Featuring a new
deck design with
more leg room!



New and improved deck design for increased comfort and leg room.



Intuitive control panel framed by a soft-touch delta tiller handle.



EZ Turn technology provides unrivaled tight turning without sacrificing stability.

FEATURES

- **New sloped deck** design for improved comfort and leg room
- **Front and rear pathway lights**
- **EZ Turn** technology provides unrivaled tight turning without sacrificing stability
- **Feather-touch disassembly** allows for easy separation into 5 lightweight pieces for transport or storage
- Exclusive **Charge360** magnetic charging port requires minimal dexterity and lights up to signal a solid connection
- Tiller angle can be quickly adjusted with a convenient foot pedal
- Large front basket included
- Interchangeable shroud panels can be replaced with accessory colors shown below

INCLUDED COLOR PANELS (MATTE)



ACCESSORY COLOR PANELS (MATTE)



ACCESSORIES*

- Additional shroud panel colors
- Cup holder
- Military patches
- Oxygen holder
- Rear basket
- Rear mount 3" double cane/crutch holder
- Rear view mirror
- Saddle bag
- Safety flag
- Single rear mount 3" cane/crutch holder
- Under seat storage
- Walker holder
- Weather cover

* Only one rear accessory can be used at a time.