

Features, Performance, and **Travel Ready**



4-wheel

GO GO
ELITE TRAVELLER® PLUS



3-wheel

The **Go-Go Elite Traveller® Plus** offers style and performance for your on-the-go lifestyle. Standard features include a delta tiller with wraparound handles, feather-touch disassembly, black, non-scuffing tires, a front basket and an 18 AH battery pack. A 300 lb. weight capacity and a maximum speed up to 4.5 mph complement the versatility of the Go-Go Elite Traveller® Plus.





Feather-Touch Disassembly

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

| pridemobility.com |



Specifications 3W: SC53 4W: SC54

| | |
|-------------------------------------|--|
| Weight Capacity | 300 lbs. |
| Maximum Speed ¹ | Up to 4.5 mph |
| Ground Clearance ² | 1.7" at transaxle, .83" at anti-tip, 2.6" at side deck |
| Turning Radius ² | 3W: 34.5" 4W: 45.75" |
| Overall Length ² | 3W: 40" 4W: 40.25" |
| Overall Width ² | 21.25" |
| Front Tires | 3W: 2" x 8" solid 4W: 3" x 9" solid |
| Rear Tires | 3" x 9" solid |
| Range Per Charge ^{1,4} | 3W: Up to 12.7 miles 4W: Up to 12.7 miles |
| Total Weight Without Batteries | 3W: 86 lbs. 4W: 87 lbs. |
| Heaviest Piece When Disassembled | 3W: 32 lbs. (rear section) 4W: 34 lbs. (front section) |
| Standard Seating | Type: Foldable, compact, molded plastic Weight: 23 lbs. Material: Black Vinyl Dimensions: Width: 17" Depth: 17" |
| Motor Output | 24-volt DC motor, sealed transaxle |
| Brakes | Regenerative and electromechanical |
| Battery Requirements ^{3,5} | Size: (2) 18 AH Weight: 14 lbs. each 29.5 lbs. per battery pack |
| Battery Charger | Off-board, 3.5 amp |
| Colors | 2 sets with every scooter!   |

Red Blue

HCPCS Code: K0806 (Group 2, Standard)



Pride offers a variety of lifts to **transport your mobility product with ease.** To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at pridemobility.com/lift-finder/

Features

- Charger port on tiller
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Delta tiller with wraparound handles
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability
- Front tiller-mounted basket
- Includes two sets of easily changeable, red and blue colored panels
- Long-lasting, high-intensity, front LED headlight offers optimal pathway illumination
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scuffing tires
- Removable footboard basket for additional storage (4-wheel version only)

Accessories

- Additional 18AH battery box
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 50" Lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage (assembly required)
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger

Option

- Additional 18AH battery box

Military Patches



- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ³ AGM or gel cell type required.
- ⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.