

<b>Part Number</b>	SC54ALX
<b>Weight Capacity</b>	325 lbs.
<b>Maximum Speed</b> <sup>2</sup>	Up to 4.8 mph
<b>Ground Clearance</b> <sup>4</sup>	1.39" at motor
<b>Turning Radius</b> <sup>4</sup>	38.39"
<b>Overall Length</b> <sup>4</sup>	42.48"
<b>Overall Width</b> <sup>4</sup>	20.5"
<b>Seat-to-Deck Height</b>	16.18" - 17.72"
<b>Seat-to-Ground Height</b>	21.85" - 23.54"
<b>Suspension</b>	Full front and rear suspension
<b>Tires</b>	Front: 2"x 8" Rear: 2.5"x 8"
<b>Range Per Charge</b> <sup>2,3</sup>	16 Ah: Up to 18 miles at 200 lbs. 16 Ah: Up to 13.6 miles at 325 lbs. 8 Ah: Up to 9 miles at 200 lbs. 8 Ah: Up to 6.8 miles at 325 lbs.
<b>Product Weight</b>	83.2 lbs. without batteries
<b>Heaviest Piece when Disassembled</b> <sup>4</sup>	35.2 lbs.
<b>Standard Seating</b>	Type: Mesh Back Material: Black Vinyl Dimensions: 17"x17" Weight: 23 lbs.
<b>Motor Output</b>	24-volt DC motor, sealed transaxle
<b>Brakes</b>	Regenerative and electromechanical
<b>Battery Requirements</b> <sup>5</sup>	Type: Lithium battery box Size: 8 Ah or 16 Ah
<b>Battery Weight</b>	8 Ah: 8 lbs. (3.5 kg) 16 Ah: 12.2 lbs. (5.5 kg)
<b>Battery Charger</b>	Off-board 3.5 amp lithium battery charger
<b>Warranty</b>	Lifetime limited on frame; 2-year limited on electronics; 2-year limited on drivetrain 2-year limited on battery

## HCPCS Code: K0899

- <sup>1</sup> Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- <sup>2</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- <sup>3</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- <sup>4</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- <sup>5</sup> Battery weight may vary based on manufacturer.



**GOGO**<sup>®</sup>  
**ENDURANCE AL+**

Lightweight mobility **plus**  
so much more



**PRIDE**  
**MOBILITY**<sup>®</sup>

800-800-8586 • [pridemobility.com](http://pridemobility.com)



Lightweight  
aluminum  
frame



Unrivaled  
tight turning



Charge360  
magnetic  
charging



Front  
and rear CTS  
suspension

\*K0899 Power Mobility Device not coded by DME PDAC. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

081425

**PRIDE**  
**MOBILITY**<sup>®</sup>



[pridemobility.com](http://pridemobility.com)

The Go Go® Endurance AL+ is an FDA Class II Medical Device designed to aid individuals with mobility impairments.

# GO GO<sup>®</sup> ENDURANCE AL+

## Lightweight mobility **plus** so much more



FDA Class II  
Medical Device<sup>1</sup>

### FEATURES

- New **ALUMINUM FRAME** for lighter weight and improved durability
- **EZ TURN** technology for the tight turning of a 3-wheel with the stability of a 4-wheel
- **CHARGE360** allows you to magnetically charge your scooter, requiring minimal dexterity and no need to line up any pins
- **INFINITE TILLER ADJUSTMENT** lever can be operated with your hand or foot
- Front and rear pathway lights help you see and be seen in dark settings
- Feather-touch disassembly allows the unit to be easily separated into 5 pieces
- 17" x 17" mesh back seat with adjustable flip-up armrests
- USB charger located in the tiller
- Large front basket (standard)
- 8 Ah airline-safe lithium-ion battery pack (16 Ah battery is available as an optional upgrade for increased range) \*\*
- Interchangeable color panels for front and back:

Garnet Red

MATTE

Euro Grey

### OPTIONAL ACCESSORIES \*

- Cell phone holder
- Cup holder
- Double cane crutch holder
- Lap belt
- Military patches
- Oxygen tank holder
- Rear basket
- Rearview mirror (single or pair)
- Saddlebag
- Safety Flag
- Single cane crutch holder
- Under seat storage
- Walker holder
- Weather cover



\* Cup holder, rear basket, and saddlebag are available as a bundle. A maximum of one rear accessory can be installed at a time. Compact, foldable seat available in 20x18 in place of 17x17 mesh seat.

\*\* Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Consult your airline prior to flight to ensure your lithium battery is permitted.

**CHARGE360**  
magnetic charging

Pride's patented  
**EZ TURN** technology

**INFINITE TILLER  
ADJUSTMENT** can be  
operated by hand or foot

Front and rear **CTS**  
suspension