

# THE ULTIMATE, COMPACT, FOLDABLE TRAVEL SCOOTER



**GO GO**®  
**FOLDING SCOOTER**




Ready for transport and easily folded, the 4-wheel **Go-Go® Folding Scooter** is perfect for the active lifestyle. Offering great performance indoors and out at a maximum speed of 4 mph and a range per charge of up to 12.5 miles, the Go-Go Folding Scooter provides a comfortable and convenient ride. It is a super portable space saver that adjusts to fit perfectly inside any closet, corner, or vehicle trunk. The Go-Go Folding Scooter is sleek, stylish, and extremely easy to operate, making it the perfect travel mobility scooter choice.

PRIDE® **MOBILITY SCOOTERS**  
LIVE YOUR BEST®

| [pridemobility.com](http://pridemobility.com) |



### Specifications 4W: S19

Weight Capacity	250 lbs.
Maximum Speed <sup>1</sup>	Up to 4 mph
Ground Clearance <sup>2</sup>	2.3"
Turning Radius <sup>2</sup>	32"
Overall Length <sup>2</sup>	37"
Overall Width <sup>2</sup>	19"
Seat-to-Ground Height	20.5"
Seat-to-Deck Height	13"
Front Tires <sup>2</sup>	1.5" x 7" solid
Rear Tires <sup>2</sup>	2" x 8" solid
Suspension	N/A
Range Per Charge <sup>1,3</sup>	Up to 13 miles ( <i>lithium-ion</i> ); Up to 9.3 miles ( <i>sealed lead acid</i> )
Total Weight Without Batteries	51 lbs.
Heaviest Piece When Disassembled	N/A (unit doesn't disassemble)
Standard Seating	Type: Foldable, compact, molded plastic Material: Gray Vinyl Dimensions: 16" x 13"
Drive System	4-Pole 270w
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>4</sup>	Size: (2) 12 AH ( <i>sealed lead acid</i> ) standard Weight: 9 lbs. each (2 required), <b>or</b> Size: (1) 11.5 AH ( <i>lithium-ion</i> ) optional Weight: 5 lbs. each (1 required)
Battery Charger	2-amp, off-board (sealed lead acid) or 2.5-amp, off-board (lithium-ion)
Colors	 White

### HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC  
OR DOES NOT MEET CRITERIA

(US) 800-800-8586

(CAN) 888-570-1113

pridemobility.com



## Features

- Unique folding design
- Folds in just three easy steps
- Height adjustable tiller
- Up to 12.5 miles range per charge
- Easy to transport at home, hotels, or airports
- Can be stowed upright or laid down
- Articulated front end adds greater stability
- Front and rear LED lights
- Optional armrests
- Folded: 30"H, 19"W, 18"D

### To Fold:



Fold Seat



Pull Up Lever to Fold



Fold and Lock Tiller

### To Unfold:



Pull Lever



Unfold Scooter



Unfold Seat; Lock in Place



**Note:** Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go-Go<sup>®</sup> Folding Scooter and its lithium-ion battery is permitted.

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 2 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 3 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.
- 4 Battery weight may vary based on manufacturer.

All specifications provided by Heartway<sup>®</sup> Medical Products.  
Product distributed by Pride Mobility<sup>®</sup>.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.