

Part Number	Jazzy® Carbon 27
Weight Capacity	300 lbs.
Maximum Speed ¹	Up to 3.7 mph
Ground Clearance	1" at motor
Turning Radius ³	33.4"
Overall Length ⁴	35.04"
Overall Width ⁴	23.43"
Suspension	Full front suspension only
Tires	Front: 6.89" Rear: 8.27"
Range Per Charge ^{1,2}	Up to 9.3 miles per charge
Product Weight ⁵	27.8 lbs. without batteries
Standard Seating ³	Dimensions: 17" x 17" seat sling with removable cushion
Motor Output	24V 250W Brushless
Brakes	Regenerative and electromechanical
Battery Requirements	Type: Lithium battery box Size: 12 Ah
Battery Weight ⁴	4.6 lbs
Battery Charger	2 Amp
Warranty	5 years on frame; 1-year limited on electronics 1-year limited on drivetrain 1-year limited on battery

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Includes standard seat weight.
- 6 Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® and its lithium-ion battery are permitted.

* Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

**PRIDE
MOBILITY®**
800-800-8586 • pridemobility.com



Jazzy®
carbon 27

Only
27 lbs.
without
the battery



FDA Class II
Medical Device*

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



Folds and
unfolds for easy
transport



Quick disconnect
joystick with
LCD display



Full Front
suspension

Jazzy
carbon 27

Only
27 lbs.
without
the battery

Shown with **Optional** headrest
and attendant control bracket

17" x 17" seat
with mesh back

Adjustable lap belt
(Standard)

Durable, lightweight
carbon fiber frame

Mesh cup holder
(Standard)

Compact base
folds & unfolds
for easy transport

FDA Class II
Medical Device*

Full front
suspension



Folds & unfolds
for easy transport



Quick disconnect joystick
with LCD display



Full front
suspension



Accessible
release lever

PUSHING THE BOUNDARIES OF LIGHTWEIGHT MOBILITY

World-class Jazzy® quality meets the next level of lightweight mobility with the Jazzy® Carbon 27 power wheelchair. Weighing 27.8 lbs. without the battery, this unit folds in a flash, is easily stored for convenient transport, and has a durable carbon fiber frame. Travel worry-free and experience the Jazzy® difference.

FEATURES

- Unit weighs 27.8 lbs. without the battery
- 300-lb. weight capacity
- Durable carbon fiber frame
- Front suspension for a smooth ride
- Top speed of up to 3.7 mph
- Folded dimensions: 27.5" L x 23.4" W x 13.38 H
- Easy fold design secures with a strap
- Flip-up lightweight carbon fiber arms
- 12 Ah, airline compliant⁶ lithium battery (standard)
- Quick disconnect LCD joystick with digital display, USB-C port & puddle light
- Mesh cup holder, cushion and adjustable lap belt (standard)

COLOR OPTIONS

Night Sky

Black Cherry

OPTIONAL ACCESSORIES

Additional 12 Ah airline-compliant lithium battery⁶
Additional 20 Ah lithium battery (COMING SOON)
Battery travel bag
Headrest and attendant control bracket
Jazzy® travel kit
Simplicity cushion – 17" x 17"



Jazzy® travel kit shown