

<b>Part Number</b>	Jazzy® Carbon 27X
<b>Weight Capacity</b>	275 lbs.
<b>Maximum Speed</b> <sup>1</sup>	Up to 3.7 mph
<b>Ground Clearance</b> <sup>3</sup>	3"
<b>Turning Radius</b> <sup>3</sup>	35.4"
<b>Overall Length</b> <sup>3</sup>	41"
<b>Overall Width</b> <sup>3</sup>	23"
<b>Seat-to-Ground Height</b>	19.7"
<b>Tires</b>	Front: 7" Rear: 10"
<b>Range Per Charge</b> <sup>1,2</sup>	Up to 8 miles per charge
<b>Product Weight</b> <sup>3,5</sup>	27.5 lbs. without batteries
<b>Standard Seating</b> <sup>3</sup>	Cushioned back with removable seat cushion
<b>Folded Dimensions</b>	34" L x 12" W x 30" H
<b>Motor Output</b>	24V 150W
<b>Brakes</b>	Regenerative and electromechanical
<b>Battery Requirements</b> <sup>6</sup>	24V 10Ah
<b>Battery Weight</b> <sup>4</sup>	3.5 lbs
<b>Battery Charger</b>	3.5A
<b>Warranty</b>	5-year limited on frame; 13-month limited on electronics 13-month limited on drivetrain 13-month on battery

### HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

4 Battery weight may vary based on battery manufacturer.

5 Includes standard seat weight.

6 Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy® and its lithium-ion battery are permitted.

\* Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.

# Jazzy® carbon 27X

Weighs just  
**27 lbs**  
without the  
battery

## PRIDE MOBILITY®

800-800-8586 • pridemobility.com



HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.



FDA Class II  
Medical Device\*



FDA Class II  
Medical Device\*

# Jazzy<sup>®</sup> carbon 27X

Weighs just  
**27lbs**  
without the  
battery



Programmable quick disconnect joystick

Airline compliant<sup>4</sup> lithium battery (10Ah only)

Durable, lightweight carbon fiber frame

Adjustable foot rests

FDA Class II Medical Device\*

Storage pocket in seat back

Adjustable lap belt (Standard)

Removable cushion (Standard)

Carbon fiber cup holder (Standard)

X-frame design folds vertically for easy transport

## YOUR MOBILITY X-FACTOR.

The Jazzy<sup>®</sup> Carbon 27X folds vertically in just a few seconds thanks to its X-frame design, and it weighs a mere 27.5 lbs. without the battery. Comprised of super-strong carbon fiber, this power wheelchair is durable enough to withstand all your adventures and light enough to keep you moving.

## FEATURES

- Unit weighs 27.5 lbs. without the battery
- 275-lb. weight capacity
- Top speed up to 3.7 mph
- X-frame design allows the unit to fold nearly flat in vertical position
- Advanced carbon fiber frame offers durability and light weight
- Folded dimensions: 34" L x 12" W x 30" H
- Flip up lightweight carbon fiber arms
- Adjustable Foot Rests
- Airline-compliant lithium battery\*\* (10 Ah only)
- Storage pocket in seat back
- Cushion and adjustable lap belt (standard)
- Programmable quick disconnect joystick (standard)
- Carbon fiber Storage with Cup Holder (standard)
- Collapsible handles for easy storage

## COLOR OPTIONS

Black



Folds & unfolds for easy transport



Advanced durable carbon fiber frame



X-frame design for vertical folding



Latch for locking unit open or closed

## OPTIONAL ACCESSORIES

- Additional 10 Ah airline-compliant battery\*\*
- Battery Travel Bag
- Jazzy Travel Kit (Coming Soon)
- Military Patches
- Simplicity Cushion 18" x 18"
- Mesh Cup Holder