

JAZZY[®] 1103 mini



Jazzy 1103 Mini shown in Candy Apple Red with gray contoured medium-back seat, foot platform and VSI controller.



THE ULTIMATE IN POWER CHAIR MANEUVERABILITY

Introducing the Jazzy[®] 1103 Mini Power Chair, providing the ultimate in indoor maneuverability. Combining the Jazzy's patented mid-wheel drive design with a shorter wheel base, the Jazzy 1103 Mini has an incredible turning radius of 18". Its compact size and sophisticated programmable electronics make it the ultimate choice for tight-space operation. The Jazzy 1103 Mini also features a fully adjustable 180° swivel seat and multiple footrest options.

Pride[®]
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • pridemobility.com

SPECIFICATIONS

Weight capacity:	250 lbs.
Maximum speed:	Up to 4 mph
Ground clearance:	3"
Turning radius:	18"
Overall length:	34.5"
Overall width:	23"
Wheels:	Drive: 10" pneumatic Casters: 8" solid Anti-tips: 6" solid
Suspension:	Limited
Drive train:	Two-motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	50 amp, PG VSI
Specialty controls:	No
Battery charger:	3 amp, on-board (standard)
Per-charge range:	Up to 25 miles
Battery requirements:	(2) 12 volt deep cycle
Battery size:	32AH, U-1
Transport weight	
Battery weight:	24.5 lbs.
Base weight:	77 lbs.
Std. seat weight:	33 lbs. (medium-back)
Seat dimensions:	20" x 18" (maximum)
Warranty:	Lifetime limited warranty on frame 2-year limited warranty on electronics 18-months limited warranty on drive motors

Colors:



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

FEATURES

- Available in all standard Jazzy colors
- Powerful two-motor, patented mid-wheel drive design offers outstanding maneuverability
- Low center of gravity design delivers a very smooth and stable ride
- 10" knobby drive wheels and optimum weight distribution provide superior traction
- Speed range adjustable up to 4 mph
- On-board battery charger

OPTIONS

- Contoured deluxe high-back seat
- Synergy™ Seat
- Solid seat pan
- Power elevating seat
- Solid tires
- Lap belt
- Cup holder
- Heavy-duty dust cover
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Rear basket
- Angle-adjustable footplates
- Saddlebag

