

# Jazzy®

## 1450

The Jazzy® 1450 offers heavy-duty construction and performance components. This bariatric power chair features a 600 lb. weight capacity, front-wheel drive for superb indoor maneuverability in the tightest of spaces and large front drive wheels for extra absorption over obstacles. The Jazzy 1450 combines superb climbing capabilities with front anti-tips for transitioning over grade changes and handling various terrains easily.



### ULTIMATE TRACTION DRIVE TIRES



**OVERSIZED REAR CASTERS  
OPTIMIZE PERFORMANCE  
AND COMFORT**



**EXTENDED SEAT SIZES  
UP TO 32" X 28"**

**AMERICA'S FAVORITE  
POWER CHAIR**

Learn more online.



[gojazzy.com](http://gojazzy.com)

## Specifications

<b>Weight Capacity</b>	600 lbs.
<b>Maximum Speed</b>	Up to 4.5 mph
<b>Ground Clearance</b>	1.5"
<b>Turning Radius</b>	26.5"
<b>Overall Length</b>	38.25"
<b>Overall Width</b>	29.25"
<b>Drive Wheels</b>	14" solid
<b>Front Wheels</b>	4" solid anti-tips
<b>Rear Wheels</b>	9" solid casters
<b>Suspension</b>	Swing Arm
<b>Drivetrain</b>	Two-motor, front-wheel drive
<b>Braking System</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard Electronics</b>	90 amp, PG VR2
<b>Battery Charger</b>	Off-board, 8 amp
<b>Per-Charge Range</b>	Up to 12 miles
<b>Battery Requirements</b>	(2) 12 volt, Group 24
<b>Battery Weight</b>	52 lbs. each
<b>Base Weight</b>	195 lbs.
<b>Standard Seat Weight</b>	52 lbs. (high-back)
<b>Standard Seat</b>	24" x 24"
<b>Max. Seat Dimensions</b>	32" x 28"
<b>Warranty (K0826/K0827)</b>	Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors
<b>Warranty (K0839)</b>	Lifetime limited on frame; 2-year limited on electronics; 18-month limited on drive motors

## Features

- Black vinyl high-back seat with height-adjustable headrest
- Flat-free drive tires
- 60" lap belt
- Heavy-duty, base mounted, flip-up height-adjustable foot platform
- Dual-layer powder coated frame for increased durability

## Options

- Synergy® Seat
- Solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Safety flag
- Heel loops
- Lap belt (70", 80", 90")
- Swing-away footrest
- Heavy-duty drop-in leg rest
- Swing-away joystick mount
- Military patches
- XLR USB charger

## Models

<b>Jazzy 1450-2VHD-SS</b>	Code: K0826
<b>Jazzy 1450-2VHD-C</b>	Code: K0827
<b>Jazzy 1450-2SPVHD-SS</b>	Code: K0839

## Colors



Candy  
Apple Red



Viper  
Blue