

Jazzy®

1450

The Jazzy® 1450 offers heavy-duty construction and performance components including upgraded motors, a reinforced frame, 90-amp PG VR2 controller and speeds up to 4.5 mph. This bariatric power chair features a 600 lb. weight capacity, front-wheel drive for superb indoor maneuverability in the tightest of spaces and large front drive wheels for extra absorption over obstacles. The Jazzy 1450 combines superb climbing capabilities with front anti-tips for transitioning over grade changes and handling various terrains easily.



*Jazzy® 1450
shown in Candy Apple Red
with high-back seat*

SPECIFICATIONS

Weight capacity	600 lbs.
Maximum speed	Up to 4.5 mph
Ground clearance	1.5"
Turning radius	26.5"
Overall length	38.25"
Overall width	29.25"
Drive wheels	14" solid
Front wheels	4" solid anti-tips
Rear wheels	9" solid casters
Suspension	Swing Arm
Drivetrain	Two-motor, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	90 amp, PG VR2
Battery charger	Off-board, 8 amp
Per-charge range	Up to 12 miles
Battery requirements	(2) 12 volt, Group 24
Battery weight	52 lbs. each
Base weight	195 lbs.
Standard seat weight	51.5 lbs. (high-back)
Standard seat	24" x 24"
Max. seat dimensions	32" x 28"
Warranty (K0826/K0827)	Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors
Warranty (K0839)	Lifetime limited on frame; 2-year limited on electronics; 18-month limited on drive motors

MODELS

Jazzy 1450-2VHD-SS	Code: K0826
Jazzy 1450-2VHD-C	Code: K0827
Jazzy 1450-2SPVHD-SS	Code: K0839

FEATURES

- Black vinyl high-back seat with height-adjustable headrest
- Flat-free drive tires
- 60" lap belt
- Heavy-duty, base mounted, flip-up height adjustable foot platform
- Dual-layer powder coated frame for increased durability

OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Safety flag
- Heel loops
- Lap belt (70", 80", 90")
- Swing-away footrest
- Heavy-duty drop-in leg rest
- Swing-away joystick mount
- Military patches
- XLR USB charger

OPTIONAL SEATING

- Synergy® Seat
- Solid Seat Pan

COLOR CHOICES



Candy
Apple Red



Viper
Blue

Pride[®]
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.