

Jazzy®

600 ES

The Jazzy® 600 ES features patented Mid-Wheel 6® technology and Active-Trac® ATX Suspension for an extraordinary level of performance, maneuverability and the tightest turning radius in the Jazzy line. The Jazzy 600 ES offers true all-around performance indoors and out making it an excellent power chair choice.



LARGE 14" DRIVE WHEELS



**SPHERICAL
OMNI-CASTERS**



**ACTIVE-TRAC®
ATX SUSPENSION**

**AMERICA'S FAVORITE
POWER CHAIR**

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	2.75"
Turning Radius	20.5"
Overall Length	35.5"
Overall Width	24.5"
Drive Wheels	14"
Front Wheels	6" solid OMNI-casters
Rear Wheels	6" solid OMNI-casters
Suspension	Active-Trac® (ATX)
Drivetrain	Two-motor, in-line, mid-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	50A, Dynamic Shark controller
Battery Charger	Off-board, 5 amp
Per-Charge Range	Up to 17 miles
Battery Requirements	(2) 12 volt, U-1
Battery Weight	25 lbs. each
Base Weight	124 lbs.
Standard Seat Weight	52 lbs. (high-back)
Standard Seat	18" x 18-20" and 20" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	Lifetime limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

Standard Color



Cardinal Red

Models

Jazzy 600 ES 2S-SS	Code: K0822
---------------------------	-------------

Jazzy 600 ES 2S-C	Code: K0823
--------------------------	-------------

*Additional charges may apply.

Features and Benefits

- Depth adjustable, limited recline high back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Seating Options*

- Deluxe contour, limited recline high back available sizes: up to 18" x 16"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 16" x 16" - 20" x 20"
- Synergy® Seat

Options

- Bell shaped footrest
- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away leg rests
- Swing-away joystick
- Unoccupied transit kit

Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger

Recommended Lifts

- Backpacker® Series
- Zeus Series
- Outlander Full Platform
- Commander 450



Find a compatible lift for your vehicle online:
PRIDEMOBILITY.COM/LIFT-FINDER

Pride
Mobility Products Corp.

(US) 800-800-8586 • (Canada) 888-570-1113 • pridemobility.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.