



Jazzy[®]

614 SERIES

The Jazzy[®] 614 Series delivers outstanding performance and superb maneuverability in a powerful, stylish package. Enhanced in-line motors, advanced ATX suspension and upgraded electronics make the Jazzy 614 Series an innovative and exceptional overall value.



Pride[®]
Mobility Products Corp.
LIVE YOUR BEST[®]

SPECIFICATIONS

Weight capacity	300 lbs. (614 HD: 450 lbs.)
Maximum speed	Up to 4.3 mph
Ground clearance	2.6"
Turning radius	22.75"
Overall length	35.75"
Overall width	25.125"
Drive wheels	14"
Front casters	6" solid
Rear casters	6" solid
Suspension	Limited
Drivetrain	Two-motor, in-line, Mid-Wheel 6 [®]
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	50 amp, controller
Specialty controls	N/A
Battery charger	Off-board, 5 amp (standard)
Per-charge range	Up to 14 miles (614 HD: up to 11 miles)
Battery requirements	(2) 12 volt, deep cycle
Battery size	NF-22
Battery weight	37 lbs. each
Base weight	121.5 lbs.
Standard seat weight	37 lbs. (medium-back) 51.5 lbs. (high-back on 614HD)
Max. seat dimensions	22" x 22" (614 HD: 24" x 24")
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

MODELS

Jazzy 614 2S-SS	(Code: K0822)
Jazzy 614 2S-C	(Code: K0823)
Jazzy 614HD 2HD-SS	(Code: K0824)
Jazzy 614HD 2HD-C	(Code: K0825)

FEATURES

- New dual enhanced in-line motors for greater efficiency, torque, range and performance
- 450 lbs. weight capacity with HD model
- 6" casters on front & rear for maximum stability
- Active-Trac[®] ATX Suspension for enhanced performance over more varied terrain
- Easy rear access to batteries
- Seating up to 24" W x 24" D with HD model
- Toolless shroud removal
- Easy serviceability
- Depth-adjustable foot platform

OPTIONS

- Solid seat pan
- Cup holder
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Synergy[®] Seat

RECOMMENDED LIFTS

Pride Lifts & Ramps

- Backpacker[®] AVP
- Backpacker[®] MV
- Backpacker[®] Plus
- Outlander
- PrideBoom 250



COLOR CHOICES



Jazzy Red



Jazzy Blue

PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[®] —

(US) 800-800-8586 • (Canada) 888-570-1113 • www.pridemobility.com