

# Jazzy®

## 614 HD

The Jazzy® 614 HD features patented Mid-Wheel 6® technology and Active-Trac® ATX Suspension for an outstanding level of performance, maneuverability and tight turning radius, all in one heavy-duty package. All-around performance indoors and out over most surfaces makes the Jazzy 614 HD an excellent heavy-duty power chair choice.



### ACTIVE-TRAC® ATX SUSPENSION



6" CASTERS



EXTENDED SEATING SIZES

AMERICA'S FAVORITE  
POWER CHAIR

Learn more online.



[gojazzy.com](http://gojazzy.com)

# Specifications

<b>Weight Capacity</b>	450 lbs.
<b>Maximum Speed</b>	Up to 4.25 mph
<b>Ground Clearance</b>	2.25"
<b>Turning Radius</b>	22.75"
<b>Overall Length</b>	35.5"
<b>Overall Width</b>	25.25"
<b>Drive Wheels</b>	14" solid
<b>Front Wheels</b>	6" solid casters
<b>Rear Wheels</b>	6" solid casters
<b>Suspension</b>	Active-Trac® (ATX)
<b>Drivetrain</b>	Two-motor, in-line, Mid-Wheel 6®
<b>Braking System</b>	Regenerative and electro-mechanical
<b>Standard Electronics</b>	50A, Dynamic Shark controller
<b>Battery Charger</b>	Off-board, 5 amp
<b>Per-Charge Range</b>	Up to 11 miles
<b>Battery Requirements</b>	(2) 12 volt, NF-22
<b>Battery Weight</b>	37 lbs. each
<b>Base Weight</b>	122 lbs.
<b>Standard Seat Weight</b>	52 lbs.
<b>Standard Seat</b>	20" x 18-20", 22" x 20"-22"
<b>Max. Seat Dimensions</b>	24" x 24"
<b>Warranty</b>	Lifetime limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

## Colors



Candy Apple Red



Viper Blue

## Models

**Jazzy 614HD 2HD-SS** Code: K0824

**Jazzy 614HD 2HD-C** Code: K0825

\*Additional charges may apply.

# Features and Benefits

- Depth adjustable, limited recline high back seat available
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Lap belt (50", 60", 70", 80", 90")

## Seating Options\*

- Deluxe contour, limited recline high back available sizes: up to 24"x24"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 16"x16" or 20"x22"
- Synergy® Seat

## Options

- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away footrests
- Swing-away joystick

## Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger

## Recommended Lifts

- Backpacker® Series
- Zeus
- Outlander Full Platform
- Commander 450



Find a compatible lift for your vehicle online:  
[PRIDEMOBILITY.COM/LIFT-FINDER](http://PRIDEMOBILITY.COM/LIFT-FINDER)



(US) 800-800-8586 • (Canada) 888-570-1113 • [pridemobility.com](http://pridemobility.com)



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.