

Jazzy®

614 HD

The Jazzy® 614 HD features patented Mid-Wheel 6® technology and Active-Trac® ATX Suspension for an outstanding level of performance, maneuverability and tight turning radius, all in one heavy-duty package. All-around performance indoors and out over most surfaces, makes the Jazzy 614 HD an excellent heavy-duty power chair choice.



ACTIVE-TRAC® ATX SUSPENSION



6" CASTERS



EXTENDED SEATING SIZES

AMERICA'S FAVORITE
POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	450 lbs.
Maximum Speed	Up to 4.3 mph
Ground Clearance	2.25"
Turning Radius	22.75"
Overall Length	35.5"
Overall Width	25.25"
Drive Wheels	14" solid
Front Wheels	6" solid casters
Rear Wheels	6" solid casters
Suspension	Active-Trac® (ATX)
Drivetrain	Two-motor, in-line, Mid-Wheel 6®
Braking System	Regenerative and electro-mechanical
Standard Electronics	50A, Dynamic Shark controller
Battery Charger	Off-board, 5 amp
Per-Charge Range	Up to 11 miles
Battery Requirements	(2) 12 volt, deep cycle, NF-22
Battery Weight	37 lbs. each
Base Weight	122 lbs.
Standard Seat Weight	52 lbs.
Standard Seat	20" x 18-20"
Max. Seat Dimensions	24" x 24"
Warranty	Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors

Colors



Candy
Apple Red



Viper
Blue

Features

- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Lap belt (50", 70", 80", 90")

Options

- Synergy® Seat
- Solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick
- XLR USB charger

Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder

Models

Jazzy 614HD 2HD-SS	Code: K0824
Jazzy 614HD 2HD-C	Code: K0825