

The future of mobility is here.

Jazzy
air[®]

EMPOWERS





The way you move through life shouldn't be limited by a power chair. **The choice is yours.**

The Jazzy Air® elevates to the primary standing height position in just 16 seconds, providing quick environmental access to countertops, table tops, cabinetry and shelving. In addition, the Jazzy Air safely drives at speeds up to 3.5 mph while elevated, allowing users to socialize with others at walking speed.

“Confidence to
keep up with my life.”

Yvette Pegues
2016 Ms. Wheelchair International



It's a...
Walk in the Park.



It's a...
Fresh Perspective.



Rise up and never miss a beat.



FEATURES...

- Achieve 10" of power adjustable seat height in 16 seconds
- Elevate or lower seat height while on the move with the "Air" button or PG VR2 joystick drive controller
- Seat and integrated footplate design provides maximum comfort in both elevated and non-elevated positions

Jazzy Air Heightens Independence.

- Shop easily
- Reach shelves
- Safely be seen in the crowd
- Enjoy concerts and sporting events
- Partake in religious services
- Enhance social interaction





It's a...
Night to
Remember.



It's...
Time Well Spent.

Personalized to be uniquely yours.

SHROUD COLORS...



Citrine
Yellow



Sapphire
Blue



Rose
Quartz



Onyx
Black



Ruby
Red



Amber
Orange



Pearl
White



Gold



Silver



Pink
Topaz

SEATING OPTIONS...

- 16" x 16"-18" depth-adjustable, high-back reclining sport seat
- Choice of 18" x 18"-20" and 20" x 20"-22" fold-flat, depth-adjustable, high-back reclining sport seat
- Choice of 16" x 16"-18", 18" x 18"-20" and 20" x 20"-22" fold-flat, depth-adjustable solid seat pan

A wide selection of gear for extra comfort.



Cell phone holder



Rear basket



Flag



Weather cover



Saddle bag



Walker holder



Oxygen holder



Cane/crutch holder



Cup holder



XLR USB charger

OPTIONS...

- Swing away joystick
- 40 Ah batteries
- Synergy cushions

COMPATIBLE LIFTS...

- Backpacker® Plus
- Backpacker® AVP
- Outlander FP
- Outlander LP



To find a compatible lift for your vehicle and mobility product, use our new tool: pridemobility.com/lift-finder/



It's...
Mobility Made Easy.



Jazzy Air Empowers.

Whether at home or on the town, the unique appeal of Jazzy Air® will elevate user's independence and social confidence in style.

US: 800.800.8586 | CAN: 888.570.1113 | pridemobility.com |    