

Jazzy[®]

ELITE 14

The Jazzy[®] Elite 14 is designed for excellent indoor maneuverability and outstanding outdoor capability. The Jazzy Elite 14 features 14" drive wheels for enhanced outdoor performance combined with a front-wheel drive design to provide tight turns around corners. Large front wheels offer extra absorption over obstacles for superb climbing capabilities along with front anti-tip wheels for transitioning over grade changes and handling various terrains with ease.



LARGE 14" DRIVE WHEELS



PERFORMANCE-MATCHED
ELECTRONICS



COLOR-THROUGH SHROUD
MEANS SCRATCHES ARE
LESS VISIBLE

AMERICA'S FAVORITE
POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	3"
Turning Radius	24"
Overall Length	35.5"
Overall Width	25.25"
Drive Wheels	14" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	50A, Dynamic Shark controller
Battery Charger	Off-board, 5 amp
Per-Charge Range	Up to 9.8 miles
Battery Requirements	(2) 12 volt, deep cycle, NF-22
Battery Weight	37 lbs. each
Base Weight	99 lbs.
Standard Seat Weight	52 lbs. (high-back)
Standard Seat	18" x 18-20" and 20" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

Colors



Jazzy Red



Jazzy Blue

Features

- Color-through shroud
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Options

- Synergy® Seat
- Solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Heel loops
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Military patches
- Swing-away joystick
- XLR USB charger

Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder

Models

Jazzy Elite 14-2S-SS

Code: K0822

Jazzy Elite 14-2S-C

Code: K0823