

Jazzy®

ELITE ES-1

The Jazzy® Elite ES-1 offers a combination of great style, comfort and performance for your active lifestyle. The in-line, front-wheel drive design provides tight turns around corners for optimal indoor maneuverability in tight spaces. Large front wheels offer extra absorption over obstacles and superb climbing capabilities combine with front anti-tip wheels for transitioning over grade changes easily. The Jazzy Elite ES-1 handles various terrains with ease.

*Jazzy® Elite ES-1
shown in Jazzy Red
with high-back seat*



Pride®
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com



SPECIFICATIONS

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	0.75"
Turning Radius	23.5"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	5" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 12 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	24.5 lbs. each
Base Weight	68.5 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

COLOR CHOICES



Jazzy
Red



Jazzy
Blue

Pride[®]
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com

MODELS

Jazzy Elite ES-11S-C | Code: K0816

FEATURES

- Color-through shroud
- Black vinyl, depth adjustable high-back seat
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick
- Height adjustable footrest bracket
- XLR USB charger

RECOMMENDED LIFTS

- Backpacker[®] Plus
- Backpacker[®] AVP
- Outlander Full Platform
- Olympian
- Commander

LIFT FINDER

Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.