


Jazzy

ELITE ES



The Jazzy® Elite ES delivers a reliable blend of power, performance and style for your active lifestyle. The in-line, front-wheel drive technology provides tight turns around corners for optimal indoor maneuverability in small spaces and large front wheels for extra absorption over obstacles.

IN-LINE FRONT-WHEEL DRIVE TECHNOLOGY



BLACK, FLAT-FREE
NON-MARKING TIRES



COMFORT HIGH-BACK SEAT
WITH HEADREST

AMERICA'S FAVORITE POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.5"
Turning Radius	24.75"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3.5 amp
Per-Charge Range	Up to 15 miles
Battery Requirements	(2) 12 volt, U-1
Battery Weight	25 lbs. each
Base Weight	69 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	16" x 16-18", 18" x 18-20", 20" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

Colors



Candy Apple Red



Viper Blue

Models

Jazzy Elite ES-2S-SS

Code: K0822

Jazzy Elite ES-2S-C

Code: K0823

*Additional charges may apply.

Features and Benefits

- Depth adjustable, limited recline comfort high back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Seating Options*

- Depth adjustable, limited recline comfort high back seat with solid seat pan, available sizes: 16"x16"-18", 18"x18"-20"
- Deluxe contour, limited recline high back available sizes: 18"x16"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 16"x16" or 18"x20"
- Synergy® Seat

Options

- Bell shaped foot platform
- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away leg rests
- Swing-away joystick

Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger

Recommended Lifts

- Backpacker® Series
- Zeus Series
- Outlander Full Platform
- Commander 450



Find a compatible lift for your vehicle online:
PRIDEMOBILITY.COM/LIFT-FINDER

Pride
Mobility Products Corp.

(US) 800-800-8586 • (Canada) 888-570-1113 • pridemobility.com •



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.