

# Jazzy®

## ELITE ES

The Jazzy® Elite ES delivers a reliable blend of power, performance and style for your active lifestyle. The in-line, front-wheel drive technology provides tight turns around corners for optimal indoor maneuverability in small spaces and large front wheels for extra absorption over obstacles. With superb climbing capabilities combined with front anti-tip wheels for transitioning over grade changes easily, the Jazzy Elite ES handles various terrains with ease.



*Jazzy® Elite ES  
shown in Jazzy Red  
with high-back seat*

## SPECIFICATIONS

<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed</b>	Up to 4 mph
<b>Ground Clearance</b>	0.75"
<b>Turning Radius</b>	23.5"
<b>Overall Length</b>	34.75" without front riggings
<b>Overall Width</b>	22.75"
<b>Drive Wheels</b>	9" solid
<b>Front Wheels</b>	5" solid anti-tips
<b>Rear Wheels</b>	6" solid casters
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, in-line, front-wheel drive
<b>Braking System</b>	Regenerative and electro-mechanical
<b>Standard Electronics</b>	40 amp, PG GC 3 controller
<b>Battery Charger</b>	Off-board, 3 amp
<b>Per-Charge Range</b>	Up to 12 miles
<b>Battery Requirements</b>	(2) 12 volt, deep cycle, U-1
<b>Battery Weight</b>	24.5 lbs. each
<b>Base Weight</b>	68.5 lbs.
<b>Standard Seat Weight</b>	43 lbs.
<b>Standard Seat</b>	18" x 18-20"
<b>Max. Seat Dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## COLOR CHOICES



Jazzy  
Red



Jazzy  
Blue

## OPTIONAL SEATING

- Solid seat pan

## MODELS

<b>Jazzy Elite ES-2S-SS</b>	Code: K0822
<b>Jazzy Elite ES-2S-C</b>	Code: K0823

## FEATURES

- Color-through shroud
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

## OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap Belt (50", 70", 80", 90")
- Elevating legrest
- Swingaway legrest
- Heel loops
- Military patches
- Height-adjustable footrest bracket
- Swing-away joystick
- XLR USB charger

## RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander
- Olympian
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

[prideliftsandramps.com](http://prideliftsandramps.com)

**Pride**<sup>®</sup>  
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • [pridemobility.com](http://pridemobility.com)



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.