

Jazzy®

ELITE ES PORTABLE

The Jazzy® Elite ES-Portable offers a combination of great style, comfort, performance and portability and disassembles into 5 convenient pieces for your active lifestyle. The in-line, front-wheel drive design allows tight turns around corners, with large front wheels for extra absorption over obstacles.



CONVENIENT DISASSEMBLY



PERFORMANCE-MATCHED ELECTRONICS



COMFORT HIGH-BACK SEAT WITH HEADREST

AMERICA'S FAVORITE
POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.5"
Turning Radius	25.5"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 11 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	25 lbs. each
Base Weight	48 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

Colors



Candy
Apple Red



Viper
Blue

Features

- Easy disassembly for convenient transport
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Options

- Solid seat pan
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Height-adjustable footrest bracket
- Lap belt (50", 70", 80", 90")
- Cup holder
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Swing-away joystick
- Military patches
- XLR USB charger

Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder

Models

Jazzy Elite ES-P-SS

Code: K0820

Jazzy Elite ES-P-C

Code: K0821