

Jazzy[®]

SELECT 6

The Jazzy Select[®] 6 incorporates sporty design and performance features for a truly unique ride experience. Pride's patented Active-Trac[®] Suspension combined with six wheels on the ground and an innovative articulating beam allow for superior performance indoors and out, with added stability over uneven terrain and ramp transitions. Mid-wheel drive provides the tightest turning radiuses in their class and caster wheels absorb the shock of obstacles, offering outstanding handling and control.



*Jazzy Select[®] 6
shown in Jazzy Red with
high-back seat and Jazzy[®]
Armor[®] (standard)*

Pride[®]
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com



SPECIFICATIONS

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.75"
Turning Radius	21.87"
Overall Length	35.375"
Overall Width	23.5"
Drive Wheels	10"
Front Wheels	6" solid casters
Rear Wheels	6" solid casters
Suspension	Active-Trac® (ATX)
Drivetrain	Two-motor, in-line, mid-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 15 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	24.5 lbs. each
Base Weight	100 lbs.
Standard Seat Weight	43 lbs. (comfort, high-back)
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

COLOR CHOICES



Jazzy
Red



Jazzy
Blue

OPTIONAL SEATING

- Synergy® Seat
- Solid seat pan

MODELS

Jazzy Select 6-SS	Code: K0822
Jazzy Select 6-C	Code: K0823

FEATURES

- Jazzy® Armor® includes shroud and controller guards to protect against daily wear and tear
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

OPTIONS & ACCESSORIES

- Power elevating seat (non-swivel)
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick
- XLR USB charger

RECOMMENDED LIFTS

- Backpacker® Plus
- Outlander Full Platform
- Backpacker AVP
- Olympian
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

prideliftsandramps.com

Pride®
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.