

Jazzy[®] 2.0¹

SELECT 6

The Jazzy Select 6[®] 2.0¹ offers a truly unique ride experience with a sporty design and superior performance. Pride's patented Active-Trac[®] Suspension combines with Mid-Wheel 6[®] drive design for a smooth, comfortable and stable ride, both indoors and out. The depth-adjustable high-back seat folds flat for greater convenience, while offering maximum comfort and support.

This foot platform offers superior adjustability for ample foot room.



FDA Class II
Medical Device¹



This convenient folding seat is depth adjustable for maximum comfort and offers replaceable back and seat covers.



AMERICA'S FAVORITE
POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed²	Up to 4.3 mph
Ground Clearance³	2.5" at footplate, 1.77" at motor bracket
Turning Radius³	23.76"
Overall Length³	43.27"
Overall Width	23.27"
Drive Wheels	10" flat free
Front Wheels	6" flat free
Rear Wheels	6" solid
Suspension	Active-Trac [®] (ATX) with Mid-Wheel 6 [®] Technology
Drivetrain	Two-motor, Mid-Wheel 6 [®] drive
Braking System	Regenerative and electromechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3.5 amp
Per-Charge Range^{3,5}	Up to 19 miles
Battery Requirements^{4,6}	(2) 12 volt, deep cycle, U-1
Battery Weight	24.5 lbs. each
Base Weight	170 lbs. without batteries
Standard Seat Weight	50 lbs. (comfort, high-back)
Weight of Heaviest Piece	114 lbs.
Max. Seat Dimensions	20" x 20" (high-back)
Seat-to-Floor Range	21.5" x 23.25"
Warranty	5-year limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

Colors



Sapphire
Blue



Ruby
Red



Pearl
White

*Additional charges may apply.

Features

- Active-Trac[®] Suspension combines with Mid-Wheel 6[®] Technology for a smooth, stable ride over uneven terrain
- Convenient, limited recline, high-back black vinyl folding seat is depth adjustable for maximum comfort and includes replaceable back and seat covers available in sizes: 16" x 16"-18", 18" x 18"-20" and 20" x 18"-20"
- Height-adjustable, flip-up foot platform offers ample foot room for maximum comfort
- Black, flat-free, non-marking tires saves your floors
- In-line joystick mount is equipped with performance-matched electronics for easy operation
- Equipped with 60" lap belt for security and peace of mind

Optional Seating*

- Limited recline, high-back black vinyl solid seat pan with depth adjustment available in sizes: 16" x 16"-18" or 18" x 18"-20"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 20" x 18" or 20" x 20"

Options

- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away footrests
- Swing-away joystick

Accessories

- Cane/crutch holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Military patches
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger

Model

Jazzy Select 6 2.0¹ | Code: K0822/K0823

- ¹ Jazzy[®] Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁴ AGM or gel cell type required.
- ⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁶ Battery weight may vary based on manufacturer.



(US) 800-800-8586 • (Canada) 888-570-1113 • pridemobility.com •



NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.