

Jazzy[®]

SPORT 2

The Jazzy[®] Sport 2 offers extra stability combined with front-wheel drive and in-line motor technology for a compact highly maneuverable package. Designed to deliver exceptional performance indoors and out.

IN-LINE FRONT-WHEEL DRIVE TECHNOLOGY



PERFORMANCE-MATCHED ELECTRONICS



COMFORT HIGH-BACK SEAT WITH HEADREST

AMERICA'S FAVORITE POWER CHAIR

Learn more online.



gojazzy.com

Specifications

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.5"
Turning Radius	24.75"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 13 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	25 lbs. each
Base Weight	69 lbs.
Standard Seat Weight	43 lbs.
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

Colors



Jazzy Red



Jazzy Blue

Features

- Jazzy® Armor® includes shroud and controller guards to protect from daily wear and tear
- Color-through shroud
- Depth adjustable, black high-back seat
- Flip-up, removable, height and width adjustable armrest
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

Options

- Cup holder
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Lap Belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Military patches
- Swing-away joystick
- Saddle bag
- Weather cover
- XLR USB charger

Recommended Lifts

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Zeus
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

pridemobility.com/Lift-Finder