

i L E V E L R E A D Y



**Q6 EDGE® HD**  
shown in *Yellow Blaze*  
with *iLevel®* technology

# Q6 edge<sup>®</sup> HD

With features like standard 4-pole motors, Mid-Wheel 6<sup>®</sup> Drive Design and ATX Suspension, the Q6 Edge<sup>®</sup> HD is engineered to meet the performance needs of the most active user.

The Q6 Edge HD accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our *iLevel®* technology.

QUANTUM<sup>®</sup>

# Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Solid Optional: 14" Pneumatic	
<b>CASTER WHEELS</b>	Front: 6" Solid Rear: 6" Solid	
<b>ANTI-TIP WHEELS</b>	N/A	
<b>MAXIMUM SPEED</b> <sup>4</sup>	Up to 4.7 mph	
<b>GROUND CLEARANCE</b>	3" (frame center)	
<b>TURNING RADIUS</b> <sup>5</sup>	19.75"	
<b>OVERALL LENGTH</b>	35.6" without foot riggings	
<b>BASE WIDTH</b>	26.25"	
<b>SEATING SIZES</b> <sup>16</sup>	<b>TB2 HD</b>	<b>TB3 HD</b>
<b>Synergy®/Static Seat</b>	<b>Width:</b> 20-28" <b>Depth:</b> 20-28"	<b>Width:</b> 20-32" <b>Depth:</b> 20-28"
<b>Power Tilt</b>	<b>Width:</b> 20-28" <b>Depth:</b> 20-28"	<b>Width:</b> 20-32" <b>Depth:</b> 20-25"
<b>Lift &amp; Tilt</b>	N/A	<b>Width:</b> 20-24" <b>Depth:</b> 20-25"
<b>SEAT-TO-FLOOR HEIGHTS</b>	<b>TB2 HD</b>	<b>TB3 HD</b>
<b>Synergy®/Static Seat</b> <sup>8</sup>	16.5-18.5"	17.5-19.5"
<b>Power Tilt</b> <sup>8</sup>	18.2-20.9"	18.5-20.5"
<b>Lift &amp; Tilt</b> <sup>8</sup>	N/A	17.75-19.75"
<b>MANUAL TILT/RECLINE</b> <sup>10</sup>	No	
<b>BATTERY SIZE</b> <sup>6</sup>	Group 24	
<b>BATTERY WEIGHT</b> <sup>14</sup>	52 lbs.	
<b>AVAILABLE ELECTRONICS</b> <sup>12</sup>	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX	
<b>BATTERY CHARGER</b>	8A Off-board	
<b>MOTOR PACKAGES</b>	4-Pole 4.7 mph <sup>1</sup>	
<b>WEIGHT CAPACITY</b> <sup>7</sup>	450 lbs. <sup>15</sup>	
<b>BASE WEIGHT</b>	143.5 lbs.	
<b>MEDICARE CODE COVERAGE</b>	Group 3: HD (K0850/K0851), Single Power HD (K0858), Multiple Power HD (K0862)	

- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) 250 lbs. weight capacity with TRU-Balance® Systems
- 10) The Synergy® Manual Tilt weight capacity is 75 lbs. and maximum dimensions are 14" W x 14" D. Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab® units can be ordered with specialty controls
- 13) Required with high-speed motors
- 14) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 15) Weight capacity is reduced to 400 lbs. with the iLevel® power adjustable seat height
- 16) A minimum seating size of 20"x20" is required with iLevel®. iLevel® maximum seating size is 24" x 25".

# Features

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® 3 HD Power Positioning
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain
- iLevel® option enables up to 10" of seat elevation while at walking speed (up to 3.5 mph)<sup>15,16</sup>

# Options

- iLevel® power adjustable seat height<sup>15,16</sup>
- TRU-Balance® 3 (450 lbs.) HD Tilt, HD Recline, and HD Tilt & Recline
- Q-Logic 2 Drive Control System
- USB mobile device charger
- Swing-away joystick mount
- LED fender lights (standard with iLevel®)



# Colors

SEAT BACKS for units without TRU-Balance® 3 are only available in *Back in Black*.



# QUANTUM®

(US) 866-800-2002 | (Canada) 888-570-1113 | quantumrehab.com |     /quantumrehab