

# PRIDE<sup>®</sup> MOBILITY

# TRAVEL POWER CHAIR<sup>1</sup>

## EXPLORE THE WORLD

Whether you are riding the rails or cruising the seas, you need a travel power wheelchair that is lightweight and portable. Consider a Pride travel power wheelchair for your next adventure.



SCAN THE QR CODE  
TO LEARN MORE ABOUT  
TRAVEL MOBILITY.



FOLLOW US ON SOCIAL MEDIA

**PRIDE<sup>®</sup>**  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

# GO CHAIR®



## FEATURES

- Easily disassembles into 5 lightweight pieces
- Dual swing-away, under seat storage bins
- Number one for portability
- Heaviest piece is only 36 lbs.
- Available in 6 colors

## ACCESSORIES

- Armrest mount cane clip
- Cane/crutch holder (3" single, 3" double)
- Cup holder
- Lap belt (50", 70")
- Oxygen tank holder
- Rear mount bracket
- Saddlebag
- Safety flag
- Weather cover
- XLR USB charger
- Additional battery pack

Ruby Red

Pearl White

Amber Orange

Citrine Yellow

Rose Quartz

Sapphire Blue

## MILITARY PATCHES



**300lbs.**  
Weight Capacity

Up To  
**3.7mph**  
Maximum Speed

**25.5"**  
Turning Radius

Up To  
**13.2**  
Miles Per Charge at 200 lbs.<sup>2</sup>

# JAZZY® PASSPORT<sup>1</sup>



## FEATURES

- Compact and lightweight design makes traveling a breeze
- Folds in just a few simple steps
- Convenient, under seat mesh storage bag and cup holder
- Airline safe lithium-ion battery available<sup>3</sup>
- Only 54 lbs. without batteries
- Available in black

## ACCESSORIES

- XLR USB charger

## OPTIONAL

- Stealth Simplicity G Cushion, 18"x18"



Black

Folded Size  
31" x 16" x 23.5"



## MILITARY PATCHES



**250lbs.**  
Weight Capacity

Up To  
**3.6mph**  
Maximum Speed

**24.25"**  
Turning Radius

Up To  
**11**  
Miles Per Charge at 200 lbs.<sup>2</sup>

## GO CHAIR® AND JAZZY® PASSPORT WARRANTY

5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 13-month limited on batteries

HCPCS Code: K0899

1. Pride® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

2. Range varies based on user weight, terrain and battery type.

3. Lithium-ion battery is subject to federal regulations and airline policy. Please consult your airline prior to flight.



**PRIDE**  
MOBILITY  
LIVE YOUR BEST®

AUTHORIZED PROVIDER

0001868-091420