

NEW L-225

FEATURES AND BENEFITS:

- Premium furniture-grade laminate/hardwood frames for added strength and durability
- 3-position full recline ideal for napping and relaxing
- Standard large dual pockets
- Easy-to-use hand control
- Removable back
- Energy-efficient transformer, features self-diagnostic electronics
- Quiet and smooth lift system
- Heavy-duty lift actuator and scissor mechanism
- Integrated battery backup in case of power failure¹

OPTIONS:

- Residential cleaning kit
- Head pillow
- Lift chair overlay
- Lumbar pillow
- Protective pad
- Non-skid leg levelers




SPECIFICATIONS	Medium
	L-225
	3-Position Chaise Lounger 5'2" - 5'10"
Weight Capacity	375 lbs.
Seat-to-Floor	19"
Seat Depth	20"
Seat Width	20"
Top of Back to Seat	26"
Back Style	Blown Fiber Split Back
Fully Padded Chaise	Yes
Heat and Massage Option	No
Seating Type	Sinuous Spring
Distance From Wall	18"
Overall Chair Width	33"
Chair Weight	111 lbs.
Warranty	2-Year Limited



AVAILABLE FABRICS

CLOUD 9




PACIFIC



STONE



WALNUT



BLACK
CHERRY

¹. Batteries not included.



Look for the Quick Ship Indicator for the shortest lead times available.

**MADE IN
USA**
OF FOREIGN & DOMESTIC COMPONENTS

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

pridemobility.com



Pride
Mobility Products Corp.
— LIVE YOUR BEST —