

AMERICA'S #1  
LIFT CHAIR BRAND

**Pride**<sup>®</sup>  
Lift Chairs

Introducing the  
**INDUSTRY'S  
FIRST**

# DUAL MOTOR SYNC SYSTEM

The Dual Motor Sync System, now standard on the LC-358XL and XXL, is an industry-first product enhancement from Pride<sup>®</sup>.

Traditional dual inline motor systems have a natural tendency to fall out of sync over time, requiring regular manual calibration. Pride's Dual Motor Sync System actively works to keep the inline motors synchronized, ensuring long lasting quality and performance.



**HERITAGE**  
COLLECTION  
LC-358XXL



To see a video demonstration of the Dual Motor Sync System, visit <http://www.pridemobility.com/liftchairs/dualmotorsync>

Note: Deluxe Heat & Massage option not available with Dual Motor Sync System. Standard Heat & Massage is available.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1-800-800-8586 (U.S.) • 1-888-570-1113 (Canada)  
[pridemobility.com](http://www.pridemobility.com) • [prideliftchairs.com](http://www.prideliftchairs.com)

**Pride**<sup>®</sup>  
Mobility Products Corp.  
— LIVE YOUR BEST —