MODEL

Oasis Collection	LC-580iM1	LC-580iL1

Position	True Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'4" - 5'9"	True Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'9" - 6'
Weight Capacity	375 lbs	375 lbs
Seat-to-Floor	20"	21"
Seat Depth	21"	21"
Seat Width	20"	21"
Top of Back to Seat	26"	29"
Back Style	Blown Fiber Biscuit Back	Blown Fiber Biscuit Back
Fully Padded Chaise	Yes	Yes
	Pocketed Coil	Pocketed Coil
Seating Type	Springs w/Memory Foam Topper	Springs w/Memory Foam Topper
Seating Type Heat & Massage		
5 71	Foam Topper	Foam Topper
Heat & Massage Distance From	Foam Topper Yes	Foam Topper Yes
Heat & Massage Distance From Wall	Foam Topper Yes 20"	Foam Topper Yes 21"

- 1 Pride® Power Lift Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- 2 Batteries not included.
- 3 Crypton® Fabric offers a five-year warranty on the moisture barrier. Does not cover punctures in the fabric.

NOTE: Measurements may vary +/- 1". Specifications subject to change without prior notice.

Recommended HCPCS: E0627–Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.



FABRIC SELECTION



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice



Oasis COLLECTION' A NEW LEVEL **OF COMFORT** FDA Class II Medical Device PWWER RECLINERS **AMERICA'S #1 POWER LIFT RECLINER** pridemobility.com l 🖸 f 🔮 🚳

WITH THE **MEMORY BUTTON HAND CONTROL**, YOUR

PERFECT POSITION

IS ALWAYS AT YOUR FINGERTIPS!



How to Set Your Memory Buttons Find your desired position Press and hold memory Button 1 and memory button 2 for three seconds. The LED lights on the memory buttons will flash to indicate that the hand control is in "set" mode.

Press and hold the desired memory button for three seconds to set the memory position.
The LED lights will stop flashing.

Repeat steps to program memory buttons 2-4. The Oasis brings a new level of comfort to Pride® Power Lift Recliners. Overstuffed padding provides elegant comfort and design, while contoured memory foam seating and pocketed coils deliver individualized support and relaxation. The Oasis provides true-infinite positioning and a standard memory button hand control.

- Premium furniture grade laminate/hardwood frames for added strength and durability
- Sinuous no-sag back springs combine durability with just the right amount of flexibility for personalized comfort
- Specially designed arms afford a luxurious feel and firm support
- Industry-exclusive quad pocket design offers extra storage – easy access top pockets provide convenient storage for your hand control and other devices
- Removable filling allows for custom lumbar support
- Standard lumbar pillow for additional support while relaxing

Options

- Deluxe heat and massage*
- Cleaning kit
- Head pillow
- Bolster pillow
- · Lithium battery backup



* Not available with memory button hand control