




LC-250 shown in
Cloud 9 Pacific and
Lexis Sta-Kleen
Chestnut fabric

Protect your investment
with the Pride®
Power Lift Recliner
Care Kit



MADE IN
USA

Classic COLLECTION

The Classic Collection from Pride® Power Lift Recliners delivers superb value and excellent comfort with features like an innovative arm design and overstuffed fiber-filled backs. In addition, Classic Collection power lift recliners feature patented, quiet and smooth operation and stylish designs to add to any home décor.

AMERICA'S #1
POWER LIFT RECLINER

POWER
LIFT
RECLINERS

US: 1-800-800-8586 • CAN: 1-888-570-1113
pridemobility.com |    

Pride
Mobility Products Corp.
— LIVE YOUR BEST! —

Classic COLLECTION

AMERICA'S #1
POWER LIFT RECLINER

POWER
LIFT
RECLINERS



LC-250 shown in
Cloud 9 Pacific fabric



SPECIFICATIONS	Medium
	LC-250
	3-Position Chaise Lounger 5'4" - 6'
Weight Capacity	375 lbs.
Seat-to-Floor	20"
Seat Depth	20"
Seat Width	20"
Top of Back to Seat	26"
Back Style	Blown Fiber Split Back
Fully Padded Chaise	Yes
Heat and Massage Option	No
Seating Type	Total Comfort
Distance From Wall	20"
Overall Chair Width	34"
Chair Weight	118 lbs.
Warranty	2-Year Limited

CLASSIC COLLECTION FABRIC SELECTION

QUICK SHIP

CLOUD 9



BLACK CHERRY



PACIFIC



STONE



WALNUT

QUICK SHIP

LEXIS STA-KLEEN® sta.kleen. PERFORMANCE FABRIC



CHESTNUT



MUSHROOM

QUICK SHIP

Look for the Quick Ship
Indicator for the shortest
lead times available.

BENEFITS

- Premium furniture grade laminate/hardwood frame for added strength and durability
- Standard head and arm covers
- Standard large dual pockets easily accessible to store items
- Removable back
- Transformer features self-diagnostic electronics with 9' cord
- Leg levelers
- Quiet and smooth lift system
- Heavy-duty lift actuator and scissor mechanism
- Fire retardant (California specifications; bulletin 117-2013)
- Integrated battery backup in case of power failure¹
- Lifetime warranty on steel lift frame
- 2-year limited warranty

OPTIONS*

- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Lift chair overlay
- Non-skid leg levelers

Additional lead time may
apply for some options.

*Additional charges apply for all options.

¹. Batteries not included.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

US: 1-800-800-8586 • CAN: 1-888-570-1113
pridemobility.com |

Pride
Mobility Products Corp.
— LIVE YOUR BEST —