



Unmatched
comfort,
performance,
and style!



USB
Charger
Port



Serenity

COLLECTION

AMERICA'S #1
LIFT CHAIR BRAND

Pride[®]
Lift Chairs

Experience the Ultimate in Lift Chair Comfort.

The Serenity Collection from Pride[®] Lift Chairs combines superior comfort, top-notch performance and unparalleled style to create a lift chair unlike any other. Featuring Pride[®] Memory Foam for enhanced cooling comfort and individually wrapped coil springs for a balanced seating experience, the Serenity Collection Lift Chair is the ultimate in lift chair comfort.

1-800-800-8586 (U.S.) • 1-888-570-1113 (Canada)
pridemobility.com • prideliftchairs.com

Pride[®]
Mobility Products Corp.
— LIVE YOUR BEST —

Serenity

COLLECTION

Experience the Ultimate in Lift Chair Comfort.

AMERICA'S #1
LIFT CHAIR BRAND

Pride
Lift Chairs



SPECIFICATIONS

	Small	Medium	Large	Petite Wide	Medium	Large
	SR-358S	SR-358M ²	SR-358L ²	SR-525PW ²	SR-525M ²	SR-525L ²
	3-Position Chaise Lounger 5'3" & Below	3-Position Chaise Lounger 5'4" - 6'0"	3-Position Chaise Lounger 5'9" - 6'2"	Infinite-Position Zero gravity Trendelenburg Chaise Lounger 5'2" & Below	Infinite-Position Zero gravity Trendelenburg Chaise Lounger 5'4" - 6'0"	Infinite-Position Zero gravity Trendelenburg Chaise Lounger 5'9" - 6'2"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	19"	20"	21"	19"	20"	21"
Seat Depth	19"	20"	22"	19"	20"	22"
Seat Width	18"	20"	22"	24"	20"	22"
Top of Back to Seat	24"	24.5"	25.5"	25.5"	26.5"	27.5"
Back Style	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Button Back	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes
Heat & Massage Option	No	No	No	No	No	No
Total Comfort Seating	Yes	Yes	Yes	No	No	No
Distance From Wall	17"	17"	17.5"	28"	29"	30"
Overall Chair Width	35"	35"	37"	38"	33"	35"
Chair Weight	118 lbs.	117 lbs.	120 lbs.	130 lbs.	126 lbs.	128 lbs.
Warranty	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated	7-Year Prorated

BENEFITS



- Industry-first hand control with USB charger port
- Standard head & arm covers
- Standard large dual pockets
- Removable back
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- Wooden frame composed of engineered furniture grade laminate/hardwood for added strength and durability
- Fire retardant (California specifications; bulletin 117-2013)
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 9' cord
- Lifetime warranty on steel lift frame
- 7-year prorated warranty

OPTIONS*


- Residential cleaning kit
- Head pillow
- Lumbar pillow
- Protective pad
- Non-skid leg levelers
- Footrest extension available (except SR-358S)
- First Up technology available on SR-358S, 358M & 358L

Additional lead time may apply for some options.

*Additional charges apply for all options.

1. Batteries not included.
2. Available with footrest extension.

DURASOFT FABRICS



Note: Measurements may vary +/- 1". Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

FIRST UP
TECHNOLOGY



Look for the Quick Ship Indicator for the shortest lead times available.

1-800-800-8586 (U.S.) • 1-888-570-1113 (Canada)
pridemobility.com • prideliftchairs.com

Pride
Mobility Products Corp.
— LIVE YOUR BEST —