

# THE MOST COMFORTABLE RECLINER IN THE WORLD!






### ATLAS PLUS 2

When you need a moment for yourself, the **Atlas Plus 2** is there for you. Customize your comfort while you watch the big game or read your favorite book. Enjoy a refreshing beverage and keep it close with the cup holder, and conveniently charge your phone with the wireless charger.

The Atlas Plus 2 has been upgraded with enhanced true-infinite positioning and an ergonomically designed, memory toggle remote that features two programmable memory position buttons, as well as up/down/recline functions and customized comfort.













## ATLAS PLUS 2



| COLLECTION          | ATLAS PLUS 21   |
|---------------------|-----------------|
| MODEL               | PLR-2985M       |
| POSITION            | True Infinite   |
| HEIGHT RANGE        | 5′4″-5′9″       |
| SIZE                | Medium          |
| WEIGHT CAPACITY     | 400 lbs.        |
| SEAT-TO-FLOOR       | 19.5″           |
| SEAT DEPTH          | 21″             |
| SEAT WIDTH          | 22"             |
| TOP OF BACK TO SEAT | 26″             |
| FULLY PADDED CHAISE | Yes             |
| FOOTREST EXTENSION  | Standard        |
| SEATING TYPE        | Total Comfort   |
| CHAIR HEIGHT        | 42"             |
| DISTANCE FROM WALL  | 16″             |
| OVERALL WIDTH       | 39″             |
| CHAIR WEIGHT        | 164 lbs.        |
| WARRANTY            | 7-Yr. Pro-rated |

1. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

\* The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.

#### INDEPENDENT FUNCTIONS



1. The power headrest lifts your head, neck and shoulders for maximum support.



2. The upward lift offers convenient stand assist for getting in and out of the recliner.



3. The power lumbar enables you to personalize your comfort.



**4.** The footrest extension allows you to relax and stretch out completely.



**5.** The backrest moves up and down independently to help you find the perfect position.



6. **NEW!** Enhanced true-infinite position recline enables the backrest and footrest to move independently so that a variety of positions can be reached including raising your feet above your heart.

### STANDARD FFATURES

- 7. NEW! Ergonomically designed toggle remote has two programmable memory position buttons. Performs the up/down, recline functions. Additionally, it has a convenient USB charger port that allows you to charge your smart device.
- 8. Lithium battery backup ensures your recliner is operational when the power goes out.

STANDARD FABRIC: BADLANDS





