

THE MOST COMFORTABLE RECLINER IN THE WORLD!





ATLAS

The **Atlas Collection** by VivaLift® delivers comfort and relaxation that is fully customizable to meet your needs. A convenient remote with soft-touch buttons adjusts each function independently. Enjoy a beverage with the accessible cup holder. The Atlas collection offers an array of elegant fabrics with contrast stitching to complement any living room or family room.



















COLLECTION	ATLAS ¹
MODEL	PLR-985M
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift 5'4"-5'9"
SIZE	Medium
WEIGHT CAPACITY	400 lbs.
SEAT-TO-FLOOR	19.5″
SEAT DEPTH	21″
SEAT WIDTH	22″
TOP OF BACK TO SEAT	26"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	16"
OVERALL WIDTH	39″
CHAIR WEIGHT	161 lbs.
WARRANTY	7-Yr. Pro-rated

1. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

INDEPENDENT FUNCTIONS



1. The power headrest lifts your head, neck and shoulders for maximum support.



2. The upward lift offers convenient stand assist for getting in and out of the recliner.



3. The power lumbar enables you to personalize your comfort.



4. The footrest extension allows you to relax and stretch out completely.



5. The backrest moves up and down independently to help you find the perfect position.



6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.



STANDARD FEATURES

- 7. The USB remote adjusts all functions independently and allows you to charge your smart device.
- 8. Lithium battery backup ensures your recliner is operational when the power goes out.

STANDARD FABRIC: BADLANDS CONTINUES STANDARD FABRICS BADLANDS





(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

