

VivaLIFT!® ...BECAUSE YOU DESERVE MORE THAN JUST ORDINARY!



COLLECTION	PREMIER ^{1,2,3}	ULTRA ^{1,2,3}	ULTRA ^{1,2,3}	ULTRA ^{1,2,3}	RADIANCE ^{1,2,3}	RADIANCE ^{1,2,3}	RADIANCE ^{1,2,3}	RADIANCE ^{1,2,3}	TRANQUIL 2 ¹	TRANQUIL 2 ¹	TRANQUIL 2 ¹	METRO 2 ¹
MODEL	PLR-5900M	PLR-4955S	PLR-4955M	PLR-4955L	PLR-3955PW	PLR-3955S	PLR-3955M	PLR-3955LT	PLR-935S	PLR-935M	PLR-935LT	PLR-925M
SIZE	Medium	Small	Medium	Large	Petite Wide	Small	Medium	Large/Tall	Small	Medium	Large/Tall	Medium
POSITION*	True-Infinite with Tilt 5'5" - 6'1"	True-Infinite with Tilt 5'4" and below	True-Infinite with Tilt 5'5" - 6'	True-Infinite with Tilt 6' and above	True-Infinite 5'4" and below	True-Infinite 5'4" and below	True-Infinite 5'5" - 6'1"	True-Infinite 6'2" and above	True-Infinite 5'3" and below	True-Infinite 5'4" - 5'9"	True-Infinite 5'10" and above	True-Infinite 5'3" and below
OVERALL WIDTH	39.5"	37"	39"	39"	41"	37"	39"	40"	33.75"	37.5"	38"	33"
CHAIR HEIGHT	43"	41"	43"	47"	40"	41"	43"	49"	40"	42"	47"	41"
SEATING TYPE	Pocketed Coiled Springs	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort
SEAT WIDTH	23"	20"	22.5"	23"	24"	20"	22.5"	23.5"	20"	22"	23.5"	21"
SEAT DEPTH	21"	18"	21"	22"	18"	18"	21"	22"	18.5"	20.5"	22.5"	21"
SEAT TO FLOOR HEIGHT	20.5"	18"	19.5"	21.5"	20.5"	18"	19.5"	21.5"	18"	19.5"	21.5"	19"
TOP OF BACK TO SEAT	26"	25"	28"	30"	26"	25"	28"	31"	25"	27"	30"	26.5"
DISTANCE FROM WALL	17"	15"	15"	21.5"	17"	15"	15"	21.5"	13"	16"	16"	14"
HEAT	Standard	Standard	Standard	Standard	Standard	Standard	Standard	Standard	N/A	N/A	N/A	N/A
STANDARD FABRICS	SORRENTO	CAPRICCIO	SORRENTO	TRENTO ⁴	CANYON				ASTRO			SAVILLE

All VivaLift!®
Models Feature:

- Power Headrest
- Charging Port
- Power Lumbar
- Lithium Battery Backup
- Power Backrest
- 7 year Pro Rated
- Foot Extension

OUR COMMITMENT TO QUALITY IS UNMATCHED.

1. FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.
3. Unit has three different heating pads in the shoulders, lumbar and seat.
4. Available for the PLR-4955M and PLR-4955L only

*Height ranges should be used as a guide only due to varying factors of each individual user. Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.



FDA Class II Medical Device¹

VivaLIFT! PREMIER

Immersive Audio Speakers
Speakers are built into this chair & connect to TV/Radio/Media via Bluetooth® tech.

Hand Controller Features
Magnetic docking, memory button and tilt control. Connection in both pockets.

Convenience Features
Hideaway wireless charger, self-leveling cup holder, USB-A and USB-C charging ports.

Shown in Sorrento Leather Eclipse

FDA Class II Medical Device¹

VivaLIFT! ULTRA

Available in 100% American Leather
Seven fabric/color options including Sorrento leather, Trento cloth, and Capriccio.

Heat and Massage
Delivering gentle pressure on your lower back, & melting away daily stress.

Exclusive Tilt Technology
Recline your chair into a position where your knees are elevated above your heart, providing the feeling of weightlessness.

Shown in Trento⁴ Brown

My VivaLIFT! APP

Control your **VivaLIFT!** from your smartphone or tablet with ease!

Download on the App Store | GET IT ON Google Play

VivaLIFT! RADIANCE

Foot Extension

Heating Pads

FDA Class II Medical Device¹

Shown in Canyon Steel

VivaLIFT! TRANQUIL 2

FDA Class II Medical Device¹

Tranquil 2 shown in Astro Mushroom

VivaLIFT! METRO 2

FDA Class II Medical Device¹

Metro 2 shown in Saville Brown

Foot Extension

Lithium Battery Backup

USB Charger in Remote

Power Headrest