

VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

NOW
AVAILABLE

TRUE-INFINITE
POSITIONING AND
MEMORY TOGGLE
HAND CONTROL!



Shown in Badlands Walnut Fabric



COUNTLESS
POSITION
COMBINATIONS

ELEGANCE 2¹

A truly elegant, sophisticated recliner is more than just physically alluring: It's the precise combination of function, style and comfort. The **Elegance 2 Collection** from VivaLift!® Power Recliners is all of these and more, providing a cozy place to relax with your family.

The Elegance 2 has been upgraded with an ergonomically designed memory toggle remote that features two programmable memory position buttons along with up/down/recline functions for customizable comfort.



Scan here for
more information
on **VivaLift!®**









PRIDE®
MOBILITY
LIVE YOUR BEST®

ELEGANCE 2¹

VivaLIFT![®]
POWER RECLINER COLLECTION¹

| COLLECTION | ELEGANCE 2 ¹ | ELEGANCE 2 ¹ |
|---------------------|-------------------------|-------------------------|
| MODEL | PLR-975M | PLR-975L |
| POSITION | True-Infinite | True-Infinite |
| HEIGHT RANGE | 5'4"-5'9" | 5'10" and above |
| SIZE | Medium | Large |
| WEIGHT CAPACITY | 400 lbs. | 400 lbs. |
| SEAT-TO-FLOOR | 20" | 20.5" |
| SEAT DEPTH | 21" | 22" |
| SEAT WIDTH | 21" | 22" |
| TOP OF BACK TO SEAT | 27" | 29" |
| FULLY PADDED CHAISE | Yes | Yes |
| FOOTREST EXTENSION | Standard | Standard |
| SEATING TYPE | Total Comfort | Total Comfort |
| CHAIR HEIGHT | 42" | 43" |
| DISTANCE FROM WALL | 14" | 16" |
| OVERALL WIDTH | 35.5" | 37.5" |
| CHAIR WEIGHT | 152 lbs. | 157 lbs. |
| WARRANTY | 7-Yr. Pro-rated | 7-Yr. Pro-rated |

INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. **NEW!** Enhanced true-infinite position recline enables the backrest and footrest to move independently so that a variety of positions can be reached including raising your feet above your heart.

STANDARD FEATURES

- NEW!** Ergonomically designed toggle remote has two programmable memory position buttons. Performs the up/down, recline functions. Additionally, it has a convenient USB charger port that allows you to charge your smart device.
- Lithium battery backup ensures your recliner is operational when the power goes out.


1. VivaLift![®] Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

STANDARD FABRIC: **BADLANDS**



MUSHROOM

WALNUT

STEEL

(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]