

THE MOST COMFORTABLE RECLINER IN THE WORLD!






ELEGANCE

A truly elegant, sophisticated recliner is more than just physically alluring: It's the precise combination of function, style and comfort. The **Elegance Collection** from VivaLift!® Power Recliners is all of these and more, providing a cozy place to relax with your family. The soft fabric and contrast stitching makes this recliner a perfect addition to your home.















ELEGANCE



| COLLECTION | ELEGANCE ¹ | ELEGANCE ¹ |
|---------------------|--|--|
| MODEL | PLR-975M | PLR-975L |
| POSITION | Infinite Lay Flat Relaxed Reclined Upward Lift 5'4"-5'9" | Infinite Lay Flat Relaxed Reclined Upward Lift 5'10" and above |
| SIZE | Medium | Large |
| WEIGHT CAPACITY | 400 lbs. | 400 lbs. |
| SEAT-TO-FLOOR | 20″ | 20.5″ |
| SEAT DEPTH | 21″ | 22″ |
| SEAT WIDTH | 21″ | 22″ |
| TOP OF BACK TO SEAT | 27" | 29″ |
| FULLY PADDED CHAISE | Yes | Yes |
| FOOTREST EXTENSION | Standard | Standard |
| SEATING TYPE | Total Comfort | Total Comfort |
| CHAIR HEIGHT | 42″ | 43" |
| DISTANCE FROM WALL | 14" | 16" |
| OVERALL WIDTH | 35.5″ | 37.5″ |
| CHAIR WEIGHT | 152 lbs. | 157 lbs. |
| WARRANTY | 7-Yr. Pro-rated | 7-Yr. Pro-rated |

1. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

INDEPENDENT FUNCTIONS



 The power headrest lifts your head, neck and shoulders for maximum support.



2. The upward lift offers convenient stand assist for getting in and out of the recliner.



3. The power lumbar enables you to personalize your comfort.



4. The footrest extension allows you to relax and stretch out completely.



5. The backrest moves up and down independently to help you find the perfect position.



6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.



STANDARD FEATURES

- 7. The USB remote adjusts all functions independently and allows you to charge your smart device.
- 8. Lithium battery backup ensures your recliner is operational when the power goes out.

STANDARD FABRIC: BADLANDS CONTINUES STANDARD FABRICS BADLANDS





MUSHROOM

WALNUT

STEEL

(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

