

VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

A TRUE-INFINITE 4-MOTOR RECLINER



FDA Class II Medical Device¹

Shown in Burbank Taupe Fabric

WITH
COUNTLESS
POSITION
COMBINATIONS

ESCAPE¹

Enjoy endless relaxation with the **Escape Collection** from VivaLift® Power Recliners. Power lumbar offers personalized comfort and a full-width power headrest lifts your head, neck and shoulders. A convenient remote with USB charger lets you power your smart device, and the industry-first Lithium battery backup ensures your lift chair functions when the power goes out. With an array of gorgeous fabrics, it's an easy decision to make this fabulous collection a part of your home.



Scan here for
more information
on **VivaLift!**®



PRIDE®
MOBILITY
LIVE YOUR BEST®







ESCAPE¹


MADE IN
USA
OF FOREIGN & DOMESTIC COMPONENTS

VivaLIFT![®]
POWER RECLINER COLLECTION¹

COLLECTION	ESCAPE ¹	ESCAPE ¹
MODEL	PLR-990iM	PLR-990iL
POSITION	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'4"-5'9"	True-Infinite Position Zero Gravity Trendelenburg Chaise Lounge 5'10" and above
SIZE	Medium	Large
WEIGHT CAPACITY	400 lbs.	400 lbs.
SEAT-TO-FLOOR	20"	20"
SEAT DEPTH	20"	22"
SEAT WIDTH	20"	23"
TOP OF BACK TO SEAT	26"	29"
FULLY PADDED CHAISE	Yes	Yes
FOOTREST EXTENSION	Standard	Standard
SEATING TYPE	Pocketed Coil Springs	Pocketed Coil Springs
CHAIR HEIGHT	44.5"	47"
DISTANCE FROM WALL	18"	21"
OVERALL WIDTH	33"	36.5"
CHAIR WEIGHT	148 lbs.	150 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated

INDEPENDENT FUNCTIONS


-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. True Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.



STANDARD FEATURES

7. The USB remote adjusts all functions independently and allows you to charge your smart device.
8. Lithium battery backup ensures your recliner is operational when the power goes out.

STANDARD FABRIC: **BURBANK**



1. VivaLift![®] Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]