

VivaLIFT!®

POWER RECLINER COLLECTION¹


THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER



FDA Class II Medical Device¹

Shown in Saville Brown Fabric



METRO 2¹

Contemporary comfort has never been more stylish! The **Metro 2** by VivaLift!® Power Recliners offers maximum relaxation with enhanced true-infinite positioning. Enjoy an ergonomically designed memory toggle remote with two programmable memory buttons and all the essentials.



Scan here for
more information
on **VivaLift!®**









PRIDE®
MOBILITY
LIVE YOUR BEST®

METRO 2¹

VivaLIFT![®]
POWER RECLINER COLLECTION¹

COLLECTION	METRO 2 ¹
MODEL	PLR-925M
POSITION	True-Infinite
HEIGHT RANGE	5'3" and below
SIZE	Medium
WEIGHT CAPACITY	400 lbs.
SEAT-TO-FLOOR	19"
SEAT DEPTH	21"
SEAT WIDTH	21"
TOP OF BACK TO SEAT	26.5"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	41"
DISTANCE FROM WALL	14"
OVERALL WIDTH	33"
CHAIR WEIGHT	142 lbs.
WARRANTY	7-Yr. Pro-rated

INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. **NEW!** Enhanced true-infinite position recline enables the backrest and footrest to move independently so that a variety of positions can be reached including raising your feet above your heart.

STANDARD FEATURES

- 7. NEW!** Ergonomically designed toggle remote has two programmable memory position buttons. Performs the up/down, recline functions. Additionally, it has a convenient USB charger port that allows you to charge your smart device.
- 8.** Lithium battery backup ensures your recliner is operational when the power goes out.

1. VivaLift![®] Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



ES = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

STANDARD FABRIC: **SAVILLE**



GREY

BROWN

(US) 800-800-8586 / (CAN) 888-570-1113 / pridemobility.com

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]