

JET™ 3

ULTRA



Jet™ 3 Ultra shown in Jet Red with contoured high-back grey vinyl seat with foot platform and PG VSI controller.

- ▶ Multiple footrest options
- ▶ Patented mid-wheel drive design provides outstanding manoeuvrability
- ▶ Speed range adjustable up to 6 km/h (4mph)
- ▶ On-board battery charger
- ▶ 25cm (10") pneumatic or flat-free drive wheels and optimum weight distribution provide superior traction
- ▶ Increased weight capacity – 136 kg. (300lbs.)


The Jet™ 3 Ultra delivers a broad range of standard performance features like Active-Trac® Suspension and a powerful 50-amp PG VSI controller on a very compact and manoeuvrable frame. Plus, the style is fresh, exciting and contemporary. These things combine with a truly economical price to make the Jet 3 Ultra an incredible overall value in personal mobility.

Features

- ▶ Compact size for superior handling in tight quarters
- ▶ Standard Active-Trac® Suspension
- ▶ 136 kg. (300 lbs.) weight capacity
- ▶ Utilizes the trapeze bar seat mount for quick adjustability and easy interchange of a broad range of seats
- ▶ Convenient front battery access
- ▶ Curved rear caster beam provides excellent stability
- ▶ Maximum speed up to 6 km/h (4mph)
- ▶ Up to 40 km/h (25 miles) of range per charge
- ▶ Standard 50-amp PG VSI controller
- ▶ 47 cm (18.5") turning radius

JET™ 3
ULTRA

Colour Options



Jet Red



Jet Blue

Specifications

Wheels:	25 cm (10") pneumatic or solid
Rear Casters:	20 cm (8") solid
Anti-tips:	12 cm (5") solid, front-mounted
Suspension:	Active-Trac®
Maximum Speed:	Up to 6 km/h (4mph)
Brakes:	Intelligent braking (electronic, regenerative disc brakes)
Ground Clearance:	6 cm (2.5")
Turning Radius:	47 cm (18.5")
Overall Size:	Length: 90 cm (35.5") Width: 58 cm (23")
Seating Options:	Contoured high-back, grey
Drive Train:	Two-motor, mid-wheel drive
Battery:	Two 12 volt, 32AH
Per-charge Range:	Up to 40 km/h(25 miles)
Battery Charger:	4 amp, on-board (standard)
Electronics:	50 amp PG VSI
Weight Capacity:	136 kg. (300 lbs.)
Weight:	Base: 44 kg. (97lbs.) Seat: 16.5kg. (36.5 lbs.)
Colours:	Jet Red and Jet Blue
Warranty:	2 year limited

Options

- ▶ Contoured seat with solid seat pan
- ▶ Lap belt
- ▶ Cup holder
- ▶ Heavy-duty dust cover
- ▶ Rear basket
- ▶ Oxygen tank holder
- ▶ Cane/crutch holder
- ▶ Walker holder
- ▶ Swing-away footrests
- ▶ Elevating leg rests
- ▶ Heel loops

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

Pride
Mobility Products
New Zealand Ltd.

38 Lansford Crescent
Avondale, Auckland 1007
New Zealand
www.pridemobility.com

Tel: 09 828 0337
Fax: 09 828 0339