



PRIDE
MOBILITY

BAJA MOBILITY

A C C E S S O R Y C O L L E C T I O N

PRIDE
MOBILITY

BAJA MOBILITY

ACCESSORY COLLECTION

Enhance your independence and mobility with our wide selection of accessories. Talk to your Authorized Pride® Provider to find out which accessories are right for you.



Weather Covers

Protect your mobility scooter from snow, rain and wind with the Pride Mobility® Scooter Weather Cover, available in three sizes. Made from durable vinyl, this weather cover is designed to protect your scooter's electronics from the elements.

SMALL: 61" (L) x 41.25" (H)

MEDIUM: 75.75" (L) x 46.5" (H)

LARGE: 84" (L) x 49" (H)



Jazzy Travel Bag Bundle*

With our Jazzy® Travel Bag bundle, your Jazzy® Carbon can be stored and transported like never before. It's designed so that you can tow your unit and batteries through airports, parking lots, and wherever else to get you from A to B. And it serves to protect your unit from damage during transportation.

* Currently only compatible with the Jazzy Carbon.



Saddle Bag

This saddle bag slides onto either armrest, featuring adjustable velcro straps and a zippered pouch on each side. Made from durable Neoprene material, the saddle bag is available in small and large sizes:

SMALL – Approx. 9.5” wide x 5.25” high x 5” deep

LARGE – Approx. 11.75” wide x 8” high x 7” deep



Front Basket

Is it difficult for you to reach items behind you? Maybe you take a lot of trips to the grocery store or pharmacy with your scooter. If this sounds like you, you may benefit from the convenience of a front basket.



Rear Basket

A rear basket for your scooter is a great option to meet your storage needs. While a saddle bag is suitable for small items, a rear basket is best for larger items. Mounted to the back of your scooter, your items are kept secure so you can drive, hands free and worry-free.



Under Seat Storage

These dual storage bins are the perfect way to store small items within arms reach. Talk to your provider to determine whether under seat storage is compatible with your scooter or power chair.



Storage Pod

Enjoy peace of mind knowing your belongings are secure throughout your journeys with our Storage Pod. Whether carrying a valuable such as a purse, or if in need of extra space for drinks, food, and other items, the Storage Pod securely shuts and locks with a key.





Baja® Mobility Cup Holder

The Baja Mobility Cup Holder easily attaches to your scooter and does all the lifting for you wherever you go, ensuring that refreshment is always within arm's reach.



Pride Mobility® Cup Holder

When you're out and about, carrying a water bottle or other container isn't always possible, which is why our cup holder is an ideal accessory to ensure an enjoyable mobility experience.



Cell Phone Holder

No matter where you go, stay connected with the RAM® X-Grip® Cell Phone Holder. When mounted to your power chair or scooter, this universal phone holder's rubber-coated tips hold your mobile device firmly in place.



Swing Away Joystick Mount

Enjoy greater maneuverability under desks, countertops and tables. Featuring a strong magnetic hold, this mount requires very little pressure to swing open. No extra mounting hardware is required for installation.



Oxygen Tank Holder

The ability to take your oxygen with you on the go keeps you happy, healthy, and mobile. Compatible with most of our scooters, this oxygen tank holder is simple to attach and easy to remove.



Lap Belt

A lap belt provides a safe and comfortable ride, especially if you encounter uneven terrain or small bumps. Our lap belts are available in several different sizes. Ask your provider which size is right for you.



XLR USB Charger

Stay connected while on the go with our XLR USB charger. This charger plugs right into the charging port of your scooter or power chair and comes with a lanyard for your convenience. Because of its charging detection feature, it switches off automatically once your device is finished charging, wasting no energy.



Safety Flag

Increase your visibility when you're out and about with our safety flag. Don't compromise on safety. This brightly colored flag is made from durable material that can handle wind and rain.



Rearview Mirrors

Our rearview mirrors make backing up safer and easier. No more straining your neck to see behind you. The mirrors attach easily to any Pride Scooter with a delta tiller. No tools needed!



Swing-away Legrests

Our swing-away legrests can be set at any angle between 60° and 90° and are fully adjustable in 5° increments. Available as a set, you can choose between two different footplate sizes: 7.25" x 6" or 8.5" x 6"



Military Patches


Proudly display your service to our country with one of these high quality embroidered patches.

Circular patches: 4" diameter
Purple Heart patch: 6" x 1.25"





Cane/Crutch Holder



Single or Double Cane Crutch Holder



Forearm Crutch Holder



Wishbone Crutch Holder



Walker Holder

Contact your provider to determine which holder is best suited for your cane, crutch, or walker and whether any mounting hardware is required.




Cane Clip Attachment

Always keep your cane within arm's reach with this convenient accessory that mounts below either armrest.



Residual Limb Supports

Please consult your provider to determine if this accessory is compatible with your power chair.



Elevated Legrests

These elevated legrests can be set at any angle between 70° and 180°. Consult your provider to determine if these are compatible with your power chair.




Heel Loops

Our heel loops are equipped with adjustable straps and will fit any of our foot rigging accessories, including our elevating legrests and swing-away legrests.



Lift Chair Head Pillow

Available in a fabric to match your lift chair, this pillow is designed to provide additional support to your head and neck.



Lift Chair Lumbar Pillow

Available in a fabric to match your lift chair, this pillow is filled with plush filling that provides firm support targeted to your lower back region.



Lift Chair Cleaning Kit

Maintain the integrity of your lift chair's fabric with this cleaning kit. Three cleaning solutions are included to remove stains along with a disinfectant cleaner that doesn't contain bleach or harsh chemicals.



Non-skid Leg Levelers

Set comes with four non-skid leg levelers. Simply screw the threaded shafts into the matching holes on the bottom of your lift chair and adjust to ensure your lift chair is on a level surface.



Lift Chair Battery Backup*

Our lithium battery backup accessory is designed to keep your lift chair operational in case of a power outage.

* All VivaLift!® models come standard with a battery backup.



Lift Chair Heat and Massage

Ask your provider if the Heat and Massage accessory is available for your lift chair.

PRIDE[®]
M O B I L I T Y


BAJA[®] M O B I L I T Y



800.800.8586



pridemobility.com



bajamobility.com

NOTE: Not all accessories are available for all products. Please speak with your Provider for details and availability. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.