

# Pride® Lifts & Ramps

Charging Innovations

## Multi-Purpose **Battery Pack**



### In-house Docking Station

- Charge your Pride Lift Battery Pack in your home by plugging into a standard 120 VAC electrical outlet
- Color changing LED shows when the Battery Pack is fully charged
- Mounting base plate makes it easy to charge anywhere

- Industry first
- Industry exclusive
- Easy installation
- Integrated smart charging system will not drain your vehicle's battery
- Charges the mobility device from the lift as you drive
- No hardwiring; just plug into auxiliary power outlet
- Trickle charge battery every time vehicle turns on
- Converts vehicle's 12v system to a mobility device's 24v system

**Pride®**  
**Lifts & Ramps**

182 Susquehanna Ave. • Exeter, PA 18643 • prideliftsandramps.com • pridemobility.com  
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

**Currently for use with Interior Products only.**