

Pride® Lifts

Mobility Made Easy.

Commander 200

Quickly load and unload your mobility device with ease using the **Commander 200**. This interior lift is the perfect user-friendly lift option for various vehicles. This lift delivers maximum maneuverability by powering up and down to load and unload your mobility device. With a 200-lb. weight capacity, this lift can be used for travel mobility, power chairs, manual wheelchairs and most foldable scooters and power chairs.

BASE WEIGHT	Trunk Mount – 43 lbs. • Boom Assembly – 16 lbs. • Mast – 9 lbs. • Base – 15 lbs.
	Hitch Mount – 49 lbs. • Boom Assembly – 16 lbs. • Mast – 12 lbs. • Hitch Mount Assembly – 21 lbs. • Hitch Tube ⁶ ◦ Class II – 8 lbs. ◦ Class III, IV – 10 lbs.
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	70"
MOTOR	12 Volt DC
WARRANTY	3-Year Transferable
HITCH CLASS	Class II, III & IV
ACCESSORY OPTIONS	Go Chair® docking kit
VEHICLE PLATFORM	Sedan, Full-size Van, Minivan, Crossover & Hatch Back
LIFTS CURRENT PRIDE® PRODUCTS	Go-Go® Travel Mobility Scooters, Jazzy® ES Series, Jazzy® Passport, Jazzy® Select Series & Jazzy® Sport 2

1. Cannot be used with hybrid vehicles.
2. Is not weatherproof, must be disassembled and stored inside of the vehicle.
3. Must be disassembled after loading / unloading mobility device.
4. Must be hard wired to vehicle battery.
5. Some disassembly of the mobility device is required.
6. Hitch mount application will come standard with one hitch tube, class II or class III / IV.

* Note: When using the trunk application the column and boom arm will need to be removed from the mount to close the trunk of the vehicle.

When using the hitch application the column and boom arm must be removed from the hitch before the vehicle is in motion. This is not an exterior lift and the hitch mount is only intended to load and unload the mobility device into the trunk of the vehicle.

To find a compatible lift for your vehicle and mobility product, use our new

LIFT FINDER

pridemobility.com/lift-finder

FEATURES

- Telescoping boom allows up to 12" of manual adjustment
- Standard single seat post docking kit or 4-point seat post karabiner kit
- Trunk mount or hitch mount available*
- Mobility seats must be removed prior to loading and unloading
- Easy assembly and disassembly in under 30 seconds
- Assemble - insert column into the base and attach the boom arm on the column
- Disassemble - remove boom arm on the column and remove column from the base
- Compatible with power chairs, scooters, manual wheelchairs, most foldable power chairs and scooters



Shown with trunk mount

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

MADE IN USA

Pride
Mobility Products Corp.
— LIVE YOUR BEST —

(800) 800-8586 (U.S.) • (888) 570-1113 (Canada) • pridemobility.com

