

# Pride® Lifts

*Mobility Made Easy*

## COMMANDER 400 Exterior Lift

The Commander 400 Exterior powered boom style lift is the perfect user-friendly lift option for pick-up trucks that provides fast loading and unloading of your mobility device. Standard with a 400 lb. lifting capacity for heavier power chairs and scooters, it is available in 2 & 3 axis. This lift delivers maximum durability, with the ability to lift and lower your mobility device to the ground over the trucks tailgate and/or over the side back into the truck bed with ease.

BASE WEIGHT*	109 lbs.
OVERALL PLATFORM LENGTH	N/A
OVERALL PRODUCT WIDTH	N/A
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	70"
MOTOR	12 Volt DC
WARRANTY	3-Year Transferable
VEHICLE PLATFORM	Full-size van, Minivan, Truck, SUV and Crossover
LIFTS CURRENT PRIDE® PRODUCTS <sup>1,2</sup>	All Pride® Scooters, Except: Pursuit® XL, Sport 36 Volt & Sport Rider  All Power Chairs Except: Jazzy® 1450 and 1650  And equivalent competitor products

### FEATURES

- 400 lb. weight capacity
- 2-axis design: Powered up/down and powered in/out
- 3-axis design: Powered up/down, powered in/out, powered adjustable telescoping boom
- Standard single seat post or 4-point seat post carabiner kit
- Coastal community powder coated
- Telescoping mounting brackets
- Complete weather-proof design
- Waterproof push-button hand control
- Compatible with power chairs and scooters



Single seat post lifting bracket

4-point seat post carabiner kit

To find a compatible lift for your vehicle and mobility product, use our new

**LIFT FINDER**

[pridemobility.com/lift-finder](http://pridemobility.com/lift-finder)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

MADE IN  
**USA**