

Pride® Lifts

Mobility Made Easy.



The **Milford Person Lift** is a highly adaptable and versatile wheelchair or scooter transfer solution. Using a sling, the MPL will transfer a person in the seated position to their vehicle, in the home, office or while traveling. MPL makes adaptations more cost effective than a 'conversion', can be easily removed and re-fitted to a new vehicle. The Milford SmartBase allows users to enjoy mobility while traveling around the home/office. The MPL can be transferred from the car to the SmartBase Unit and creates the perfect mobility lifting solution that is easy to transport, packs away in seconds, has adjustable legs and can travel anywhere, anytime.

BASE WEIGHT ¹	18 lbs.
WEIGHT CAPACITY	330 lbs. max.
MAXIMUM TRAVEL DISTANCE FROM THE GROUND	22.5"
MOTOR	12 Volt DC
WARRANTY	1-Year Transferable
HITCH CLASS	N/A
OPTIONS	<ul style="list-style-type: none">• Sling size and type• Wall bracket• Smartbase
VEHICLE PLATFORM	Full-size van, Minivan, Sedan, SUV and Crossover
LIFTS CURRENT PRIDE® PRODUCTS	N/A

1. Does not include any accessories or hardware.

FEATURES

- 330 lbs. weight capacity
- Can be fitted to the driver's side, passenger's side, front, middle or back seats
- Standard sling sizes from extra small to extra-large
- Optional amputee sling, fleece sling (for sensitive skin), net sling (for wet environments) and head support sling
- Lifts quickly and easily with the touch of a button
- Easily removed when not in use
- Vehicle applications are hardwired to vehicle
- Emergency switch and optional manual operation

MOUNTING OPTIONS

- Sill bracket & low dash adapter - use over large door sills or extended dashboards
- Dropped extension arms- available in 2",4",6",10" and 12" left and right hand drop extensions
- Wall bracket spacers- available in 3",6",9",12" extension lengths
- Spreader bar covers- soft spreader bar covers to protect the user against involuntary movements



Milford SmartBase

Optional SmartBase with removable, rechargeable battery pack

Truly adjustable transfer solutions. Quick & easy to operate

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.