

CELEBRITY X

THE NEXT GENERATION 10" CLASS SCOOTER

CELEBRITY® X


The **Celebrity® X Scooter** is a feature-rich luxury product at an attractive price. Available in a 3 or 4-wheel format, it offers elegant simplicity, unmatched reliability, impressive performance and easy-to-use operation. The Celebrity X is a perfect blend of style, performance technology, convenient standard features and affordable pricing, making it the best overall luxury scooter value in its class.



Celebrity® X 3-wheel Celebrity® X 4-wheel

Model number:	SC4001	SC4401
Weight capacity:	350 lbs.	350 lbs.
Maximum speed:	Up to 6 mph	Up to 6 mph
Ground clearance:	2.25"	2.25"
Turning radius:	41.5"	54"
Overall length:	46"	46.75"
Overall width:	24"	24"
Wheels		
Front:	10" Solid	10" Solid
Rear:	10" Solid	10" Solid
Anti-tips:	Rear	Rear
Suspension:	Rear	Front & rear
Transport weight		
Total w/o batteries:	141.5 lbs.	151 lbs.
Heaviest piece when disassembled:	56 lbs. (front section)	65.5 lbs. (front section)
Drivetrain:	Sealed transaxle, 24 volt DC motor	
Dual braking system:	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)	
Power elevating seat:	Optional: Special order required	
Std. seating:	Deluxe molded plastic with sliders	
Battery charger:	Off-board, 5 amp (standard)	
Per-charge range:	Up to 17 miles	
Battery requirements:	(2) 12V deep cycle (AGM or gel-cell type recommended)	
Battery size & weight ¹ :	40 AH: 32.5 lbs. each; U-1: 24.5 lbs. each	
Warranty:	3-year limited warranty	

Colors:



¹Battery weight may vary based on manufacturer.

FEATURES


- Delta tiller with wraparound handles
- Ultra-simple, auto-latching lockup for quick disassembly/reassembly
- Bright, long-lasting LED lights for safety
- Speed adjustable up to 6 mph
- 350 lbs. weight capacity
- Easily accessible tie down points (for transport of unoccupied scooter)
- Front frame-mounted seat post for optimal stability
- Front basket and rearview mirror standard
- Tiller rotation lock permits tiller to be used as a handle during transport

OPTIONS

- CRS deluxe reclining high-back with headrest and sliders
- Forearm crutch holder
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Saddlebag
- Quad cane holder
- Jumbo front basket
- Weather cover

Pride™ Lifts & Ramps

- Backpacker® Plus
- Outlander
- Silver Boom 300
- Backpacker® MV



These lifts are compatible with the Celebrity® X.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

PRIDE
PERFORMANCE MOBILITY
— LIVE YOUR BEST™ —