



The exciting Pride Celebrity is fully loaded with innovative standard features like Pride's exclusive, infinitely adjustable tiller with auto-leveling console that features easy to operate controls. In addition, the Celebrity offers a "soft-ride" full suspension system, a height-adjustable seat post, and full lighting package. With these features and more, the Celebrity provides unmatched convenience and performance, setting the benchmark for what a 21st Century scooter should be.

**The Celebrity  
is as stylish  
as it is high  
performing.**



**800-800-8586 (U.S.)  
888-570-1113 (Canada)  
www.pridemobility.com**

**Pride**  
Mobility Products Corp.

## Specifications

Model number(s)	3-wheel: SC4000 4-wheel: SC4400
Weight capacity	350 lbs.
Maximum speed	Up to 5.25 mph
Ground clearance	3.5"
Turning radius	3-wheel: 48" 4-wheel: 59"
Overall length	3-wheel: 50" 4-wheel: 51.75"
Overall width	24"

## Tires

Front	3" x 10" pneumatic
Rear	4" x 10" pneumatic
Anti-Tips	Rear
Standard seating	Deluxe contoured medium-back seat with sliders and pneumatic height adjustment
Suspension	Front and rear
Drive train	Sealed transaxle, 24-volt DC motor
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	On-board, 3 amp
Per-charge range	Up to 30 miles
Battery requirements	(2) 12V deep-cycle (AGM or gel cell)
Compatible battery size(s)	40AH or 32 AH, U-1

## Weight

Batteries (2 required)	40AH: 32.5 lbs. ea. 32AH: 24.5 lbs. ea.
Total w/o batteries	3-wheel: 147.5 lbs. 4-wheel: 163.5 lbs.
Heaviest piece	64 lbs. (rear)

## Colors



Candy  
Apple  
Red



Viper  
Blue

## Benefits

- Full front and rear "soft ride" suspension for a smooth and comfortable ride
- Front basket, rearview mirror and full lighting package
- All controls are at your fingertips
- On-board battery charger
- Rear accessory bracket for easy installation of rear basket, oxygen holder, cane/crutch holder and more
- Disassembles easily for transport
- External freewheel lever
- Battery/power indicator
- Easily accessible circuit breaker/charger port
- Tiller infinitely adjusts to your most comfortable driving position, while the console, mirror and headlights all remain level

## Options

- Power elevating seat
- Deluxe contoured high-back seat
- Forearm crutch holder
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- External battery charger
- Saddlebag
- Lap belt
- Cup holder
- Weather breaker canopy
- Quad cane holder
- Flat-free tires

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.