

When it comes to your busy, on-the-go lifestyle, you need a scooter that is up to the task. Pride's GO GO[®] Travel Mobility series offers a wide range of mobility scooters to suit your needs. Enjoy total convenience and exceptional features with Pride's travel scooter lineup!

Most models feature...

- Quick and easy feather-touch disassembly
- Delta tiller with port to conveniently charge your scooter
- Standard front basket

Go Go Elite
Traveller[®] 2
featuring



FDA Class II Medical Devices¹











SCAN THE QR CODE
TO LEARN MORE ABOUT
TRAVEL MOBILITY!



FOLLOW US ON SOCIAL MEDIA

PRIDE[®]
MOBILITY

		PART #	COLORS	WEIGHT CAPACITY	TURNING RADIUS ³	MAX RANGE PER CHARGE ^{2,6}	TOP SPEED ⁶	CODING ⁹
	GO GO® CARBON¹ • • • •	SC25M (4W)	Black (Carbon Fiber)	300 lbs.	53.15	Up to 10 mph	3.7 mph	K0899 Power Mobility Device, Not Coded by DME PDAC
	GO GO® SUPER PORTABLE* • • • •	SC15 (4W)	Mint, Black	300 lbs.	51.6"	Up to 10 miles at 220 lbs. Up to 8.8 miles at 300 lbs.	4.5 mph	Non-Medical*
	GO GO® ENDURANCE Li¹ • • • •	S54LXLIT	Garnet Red (matte), Euro Grey	325 lbs.	48.25"	8Ah: Up to 7.2 miles ⁴ 16Ah MAXX: Up to 14.4 miles ⁸	4.2 mph	K0899 Power Mobility Device, Not Coded by DME PDAC
	GO GO® ULTRA X¹ • • • •	S39 (3W) S49 (4W)	Black	300 lbs.	36.5" (3W) 52.25" (4W)	3W: Up to 8 miles (200 lbs.) 3W: Up to 6.2 miles (300 lbs.) 4W: Up to 8.9 miles (200 lbs.) 4W: Up to 6.6 miles (300 lbs.)	3W: 4.3 mph 4W: 4.1 mph	K0800 Group 1 Standard
	GO GO ELITE TRAVELLER®¹ • • • •	SC40E (3W)	Red, Blue	300 lbs.	33" (3W)	3W: Up to 7.25 miles (12A) 3W: Up to 10.5 miles (18A)	4 mph	K0800 Group 1 Standard
	 GO GO ELITE TRAVELLER®² • • • •	SC442E (4W)	Chell-o, White, Iceberg Blue, Sugar Plum, Robin's Egg Blue, Tangerine (all matte finish)	300 lbs.	37"	Up to 8.8 miles (12A) Up to 13.8 miles (18A)	4.23 mph	K0800 Group 1 Standard
	GO GO® SPORT¹ • • • •	S73 (3W) S74 (4W)	Red, Blue	325 lbs.	36" (3W) 50.5" (4W)	3W: Up to 12.75 miles 4W: Up to 12 miles	4.7 mph	K0801 Group 1 Heavy Duty

WARRANTIES

FRAME: Lifetime limited, 3-year limited for Go Go Super Portable
ELECTRONICS: 2-year limited; 1-year limited for Go Go Carbon
DRIVETRAIN: 2-year limited; 1-year limited for Go Go Carbon
GO GO® SLA (Sealed Lead Acid) BATTERIES: 6-month limited

GO GO® ENDURANCE Li BATTERIES: 2-year limited
GO GO® CARBON BATTERIES: 2-year limited
GO GO® SUPER PORTABLE BATTERIES: 2-year limited

ACCESSORIES⁵

- Cane clip attachment
- Cane/crutch holder
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 50" Lap belt
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger

- 1) Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2) Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3) Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4) Lithium-ion battery is subject to federal regulations and airline policy. Consult airline prior to flight.
- 5) Not all accessories are available for all models.

*This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.

- 6) Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This can be subject to a variance of (+ or -) 10%.
- 7) Red and black interchangeable shroud panels on the GO GO® Elite Traveller 2 are standard. Chell-o, White, Iceberg Blue, Sugar Plum, Robin's Egg Blue and Tangerine shroud panels are optional.
- 8) 16AH MAXX battery pack is not airline compliant.
- 9) DME suppliers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply.