



When it comes to your busy, on-the-go lifestyle, you need a scooter that is up to the task. Pride's Go-Go[®] Travel Mobility series¹ offers a wide range of mobility scooters to suit your needs. Enjoy total convenience and exceptional features with Pride's travel scooter lineup!

Most models feature...

- Quick and easy feather-touch disassembly
- Delta tiller with port to conveniently charge your scooter
- Standard front basket



FDA Class II Medical Device¹









SCAN THE QR CODE
TO LEARN MORE ABOUT
GO-GO[®] TRAVEL!¹



FOLLOW US ON SOCIAL MEDIA

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

		PART #	COLORS	WEIGHT CAPACITY	TURNING RADIUS ³	MAX RANGE PER CHARGE ^{2,6}	TOP SPEED ⁶	CODING
	GO-GO® ES 2¹	S81 (3W)	Black	250 lbs.	31"	Up to 6.8 miles	Up to 4 mph	K0800 Group 1 Standard
	GO-GO® ULTRA X¹	SC40X (3W) SC44X (4W)	Red	260 lbs.	33" (3W) 44" (4W)	3W: Up to 6.9 miles 4W: Up to 7.2 miles	Up to 4 mph	K0899 Power Mobility Device, Not Coded by DME PDAC or does not meet criteria
	GO-GO ELITE TRAVELLER®¹	SC40E (3W) SC44E (4W)	Red, Blue	300 lbs.	33" (3W) 44" (4W)	3W: Up to 8.5 miles (12A) 3W: Up to 12.4 miles (18A) 4W: Up to 8.2 miles (12A) 4W: Up to 12 miles (18A)	Up to 4 mph	K0800 Group 1 Standard
	GO-GO ELITE TRAVELLER® PLUS¹	SC53 (3W) SC54 (4W)	Red, Blue	300 lbs.	34.5" (3W) 45.75" (4W)	Up to 14.5 miles	Up to 4.5 mph	K0806 Group 2 Standard
	GO-GO® LX WITH CTS¹	SC50LX (3W) SC50LX (4W)	Red, Blue, Orion Yellow, Onyx Black, Pearl Pink, and Sunburst Orange ⁷	300 lbs.	33" (3W) 44.75" (4W)	3W: Up to 8 miles 4W: Up to 8.6 miles	Up to 4.4 mph	Power Mobility Device, Not Coded by DME PDAC or does not meet criteria
	GO-GO® SPORT¹	S73 (3W) S74 (4W)	Red, Blue	325 lbs.	36" (3W) 50.5" (4W)	3W: Up to 15.9 miles 4W: Up to 14.5 miles	Up to 4.7 mph	K0801 Group 1 Heavy Duty
	GO-GO® FOLDING SCOOTER^{1,4}	SLA: S19WH1001 Lithium-Ion S19WH1002	White	250 lbs.	32"	SLA: Up to 9.3 miles Lith-Ion: Up to 13 miles	Up to 4 mph	K0899 Power Mobility Device, Not Coded by DME PDAC or does not meet criteria

WARRANTY

GO-GO® FOLDING SCOOTER

FRAME: 2-year CONTROLLER: 18 months
ELECTRONICS & CHARGER: 1-year

ALL OTHER GO-GO® MODELS

FRAME: Lifetime limited ELECTRONICS: 2-year limited
DRIVETRAIN: 2-year limited

ACCESSORIES⁵

- Cane clip attachment
- Cane/crutch holder
- Cup holder
- Forearm crutch holder
- 50" Lap belt
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Under seat storage
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger

1) This is an FDA Class II Medical Device designed to aid individuals with mobility impairments.
2) Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
3) Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
4) Lithium-ion battery is subject to federal regulations and airline policy. Consult airline prior to flight.
5) Not all accessories are available for all models.
6) Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This can be subject to a variance of (+ or -) 10%.
7) Red and blue interchangeable shroud panels on the Go-Go® LX with CTS are standard. Orion Yellow, Onyx Black, Pearl Pink, and Sunburst Orange shroud panels are optional.

AUTHORIZED PROVIDER


pridemobility.com

PRIDE
MOBILITY
LIVE YOUR BEST®